





Features


Benefits

1	Superior front horizontal rail.	Additional rigidity. Impurity prevention.
2	Removable front center post.	Full flexibility for all door and drawer arrangements.
3	The front posts can accommodate varying arrangements of doors & drawers.	Flexibility even after installation
4	Wrap-around cabinet shell construction.	Provides maximum strength and eliminates need for a fixed finished back.
5	Removable back panel.	Allows accessibility to mechanical services.
6	Slotted rear post.	Support for slide arms and shelves. Full flexibility for future needs.
7	Cabinet body is fully enclosed.	Facilitating maintenance and cleaning.
8	Extended drawer depth.	Maximized storage capacity.
9	Insulated facades.	The doors and drawer fronts are made of ¾" thick clear maple wood.
10	Self-closing ball bearing slide arms.	Avoiding accidental spills.
11	Threaded 5/16" leveler with large foot.	Solid, sturdy and prevents scratching.
12	Various standard handles.	Adhere to design intent.

Forte Flush Overlay Door/Drawer Combination Possibilities



Our cabinet design is engineered for interchangeability of the same drawer and door heights for all different cabinet types including the following drawer and door heights:



* The top drawer and full height door on an ADA cabinet are of different heights than for all other cabinets and will only fit at this position.

Reading the product identification numbers code

The Bedcolab product identification number is simple, making it easy to verify that you have ordered the correct cabinet. Each number or letter in the designation represents an aspect of the cabinet. Below are examples of the product identification numbers you will see in this catalogue along with an explanation of each designation:

Product Identification Code





FO	A	6A	30	5
Forte Flush Overlay Identification				
Cabinet Height		Quantity of Drawers		
- A – standing - B – sitting - ADA – height		1 – 16		
		Cabinet Width		
Cabinet Style				
- Drawer cabinet - Door cabinet - Combination door/drawer - Mixed combination door/drawer		15 in. 30 in.		
		18 in. 36 in.		
		21 in. 42 in.		
		24 in. 48 in.		

Base cabinets

FO	75	42	18
Forte Flush Overlay Identification			
Cabinet Style		Cabinet Depth	
- Single door - Double door - Glass door - Glazed door		14, 18, 22* (Tall cabinets)	
		Cabinet Width	
		15 in. 30 in.	
		18 in. 36 in.	
		21 in. 42 in.	
		24 in. 48 in.	


Wall and Tall cabinets

Please refer to pages 26 and 27 of this Forte Flush Overlay catalog for hardware options and optional features.




Typical Side Views

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FO-A1-12-18	12	305
FO-A1-12-21	12	305
FO-A1-15	15	381
FO-A1-18	18	457
FO-A1-21	21	533
FO-A1-24	24	610

1 door unit




Model no.	Width	
	inches	mm
FO-A5-30-2	30	762
FO-A5-36-2	36	914
FO-A5-42-2	42	1067
FO-A5-48-2	48	1219

2 drawer & 2 door unit




FO-A2-30	30	762
FO-A2-36	36	914
FO-A2-42	42	1067
FO-A2-48	48	1219

2 door unit




FO-A5-30-3	30	762
FO-A5-36-3	36	914
FO-A5-42-3	42	1067
FO-A5-48-3	48	1219

2 half drawer, 1 drawer & 2 door unit




FO-A4-15-1	15	381
FO-A4-18-1	18	457
FO-A4-21-1	21	533
FO-A4-24-1	24	610

1 drawer - 1 door unit




FO-A5-30-4	30	762
FO-A5-36-4	36	914
FO-A5-42-4	42	1067
FO-A5-48-4	48	1219

4 half drawer & 2 door unit




FO-A4-15-2	15	381
FO-A4-18-2	18	457
FO-A4-21-2	21	533
FO-A4-24-2	24	610

2 half drawer & 1 door unit




FO-A6-30-4	30	762
FO-A6-36-4	36	914
FO-A6-42-4	42	1067
FO-A6-48-4	48	1219

4 drawer & 1 door unit




FO-A5-30-1	30	762
FO-A5-36-1	36	914
FO-A5-42-1 *3	42	1067
FO-A5-48-1 *3	48	1219

1 drawer & 2 door unit



FO-A6-30-5	30	762
FO-A6-36-5	36	914
FO-A6-42-5	42	1067
FO-A6-48-5	48	1219


2 half drawer, 3 drawer & 1 door unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.




**6 half drawer, 1 drawer
& 1 door unit**

Model no.	Width	
	inches	mm
FO-A6-30-7	30	762
FO-A6-36-7	36	914
FO-A6-42-7	42	1067
FO-A6-48-7	48	1219




**2 half drawer
& 3 drawer unit**

Model no.	Width	
	inches	mm
FO-A7-15-5	15	381
FO-A7-18-5	18	457
FO-A7-21-5	21	533
FO-A7-24-5	24	610
FO-A7-30-5 *1	30	762
FO-A7-36-5 *1	36	914
FO-A7-42-5 *1*3	42	1067
FO-A7-48-5 *1*3	48	1219




5 drawer & 1 door unit

FO-A6A-30-5	30	762
FO-A6A-36-5	36	914
FO-A6A-42-5	42	1067
FO-A6A-48-5	48	1219




**4 half drawer
& 2 drawer unit**

FO-A7-15-6	15	381
FO-A7-18-6	18	457
FO-A7-21-6	21	533
FO-A7-24-6	24	610
FO-A7-30-6 *1	30	762
FO-A7-36-6 *1	36	914
FO-A7-42-6 *1*3	42	1067
FO-A7-48-6 *1*3	48	1219




2 half drawer, 4 drawer & 1 door unit

FO-A6A-30-6	30	762
FO-A6A-36-6	36	914
FO-A6A-42-6	42	1067
FO-A6A-48-6	48	1219




**6 half drawer
& 1 drawer unit**

FO-A7-15-7	15	381
FO-A7-18-7	18	457
FO-A7-21-7	21	533
FO-A7-24-7	24	610
FO-A7-30-7 *1	30	762
FO-A7-36-7 *1	36	914
FO-A7-42-7 *1*3	42	1067
FO-A7-48-7 *1*3	48	1219




4 half drawer, 3 drawer & 1 door unit

FO-A6A-30-7	30	762
FO-A6A-36-7	36	914
FO-A6A-42-7	42	1067
FO-A6A-48-7	48	1219




8 half drawer unit

FO-A7-15-8	15	381
FO-A7-18-8	18	457
FO-A7-21-8	21	533
FO-A7-24-8	24	610
FO-A7-30-8 *1	30	762
FO-A7-36-8 *1	36	914
FO-A7-42-8 *1*3	42	1067
FO-A7-48-8 *1*3	48	1219



4 drawer unit

FO-A7-15-4	15	381
FO-A7-18-4	18	457
FO-A7-21-4	21	533
FO-A7-24-4	24	610
FO-A7-30-4 *1	30	762
FO-A7-36-4 *1	36	914
FO-A7-42-4 *1*3	42	1067
FO-A7-48-4 *1*3	48	1219




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FO-A8-30-5 *1	30	762
FO-A8-36-5 *1	36	914
FO-A8-42-5 *1*3	42	1067
FO-A8-48-5 *1*3	48	1219

5 drawer unit




Model no.	Width	
	inches	mm
FO-A10-30-10	30	762
FO-A10-36-10	36	914
FO-A10-42-10	42	1067
FO-A10-48-10	48	1219

4 half drawer & 6 drawer unit




FO-A8-30-7 *1	30	762
FO-A8-36-7 *1	36	914
FO-A8-42-7 *1*3	42	1067
FO-A8-48-7 *1*3	48	1219

4 half drawer & 3 drawer unit




FO-A10-30-12	30	762
FO-A10-36-12	36	914
FO-A10-42-12	42	1067
FO-A10-48-12	48	1219

8 half drawer & 4 drawer unit




FO-A9-30-6 *1	30	762
FO-A9-36-6 *1	36	914
FO-A9-42-6 *1*3	42	1067
FO-A9-48-6 *1*3	48	1219

6 drawer unit




FO-A10-30-14	30	762
FO-A10-36-14	36	914
FO-A10-42-14	42	1067
FO-A10-48-14	48	1219

12 half drawer & 2 drawer unit




FO-A9-30-10 *1	30	762
FO-A9-36-10 *1	36	914
FO-A9-42-10 *1*3	42	1067
FO-A9-48-10 *1*3	48	1219

8 half drawer & 2 drawer unit




FO-A10-30-16	30	762
FO-A10-36-16	36	914
FO-A10-42-16	42	1067
FO-A10-48-16	48	1219

16 half drawer unit




FO-A10-30-8	30	762
FO-A10-36-8	36	914
FO-A10-42-8	42	1067
FO-A10-48-8	48	1219

8 drawer unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FO-A11-15-2 *2	15	381
FO-A11-18-2 *2	18	457
FO-A11-21-2 *2	21	533
FO-A11-24-2 *2	24	610

2 file drawer unit




Model no.	Width	
	inches	mm
FO-A11-30-6 *2	30	762
FO-A11-36-6 *2	36	914
FO-A11-42-6 *2	42	1067
FO-A11-48-6 *2	48	1219

4 drawer & 2 file drawer unit




FO-A11-15-3 *2	15	381
FO-A11-18-3 *2	18	457
FO-A11-21-3 *2	21	533
FO-A11-24-3 *2	24	610

2 drawer & 1 file drawer unit




FO-A11-30-8 *2	30	762
FO-A11-36-8 *2	36	914
FO-A11-42-8 *2	42	1067
FO-A11-48-8 *2	48	1219

4 half drawer, 2 drawer & 2 file drawer unit




FO-A11-15-4 *2	15	381
FO-A11-18-4 *2	18	457
FO-A11-21-4 *2	21	533
FO-A11-24-4 *2	24	610

2 half drawer, 1 drawer & 1 file drawer unit




FO-A11-30-10 *2	30	762
FO-A11-36-10 *2	36	914
FO-A11-42-10 *2	42	1067
FO-A11-48-10 *2	48	1219

8 half drawer & 2 file drawer unit




FO-A11-15-5 *2	15	381
FO-A11-18-5 *2	18	457
FO-A11-21-5 *2	21	533
FO-A11-24-5 *2	24	610

4 half drawer & 1 file drawer unit




FO-A12-30-3 *2	30	762
FO-A12-36-3 *2	36	914
FO-A12-42-3 *2	42	1067
FO-A12-48-3 *2	48	1219

2 drawer, 1 file drawer & 1 door unit




FO-A11-30-4 *2	30	762
FO-A11-36-4 *2	36	914
FO-A11-42-4 *2	42	1067
FO-A11-48-4 *2	48	1219

4 file drawer unit



FO-A12-30-4 *2	30	762
FO-A12-36-4 *2	36	914
FO-A12-42-4 *2	42	1067
FO-A12-48-4 *2	48	1219

2 half drawer, 1 drawer, 1 file drawer & 1 door unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FO-A12-30-5 *2	30	762
FO-A12-36-5 *2	36	914
FO-A12-42-5 *2	42	1067
FO-A12-48-5 *2	48	1219

4 half drawer, 1 drawer, 1 file drawer & 1 door unit




FO-A13-30-4 *2	30	762
FO-A13-36-4 *2	36	914
FO-A13-42-4 *2	42	1067
FO-A13-48-4 *2	48	1219

3 drawer, 1 file drawer & 1 door unit




FO-A13-30-5 *2	30	762
FO-A13-36-5 *2	36	914
FO-A13-42-5 *2	42	1067
FO-A13-48-5 *2	48	1219

2 half drawer, 2 drawer, 1 file drawer & 1 door unit




FO-A13-30-6 *2	30	762
FO-A13-36-6 *2	36	914
FO-A13-42-6 *2	42	1067
FO-A13-48-6 *2	48	1219

4 half drawer, 1 drawer, 1 file drawer & 1 door unit




FO-A13-30-7 *2	30	762
FO-A13-36-7 *2	36	914
FO-A13-42-7 *2	42	1067
FO-A13-48-7 *2	48	1219

6 half drawer, 1 file drawer & 1 door unit




Model no.	Width	
	inches	mm
FO-A20-30	30	762
FO-A20-36	36	914
FO-A20-42	42	1067
FO-A20-48	48	1219

Corner unit



FO-A21-30	30	762
FO-A21-36	36	914
FO-A21-42	42	1067
FO-A21-48	48	1219


Corner unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.



1 door sink unit

Model no.	Width	
	inches	mm
FO-A15-15 *4	15	381
FO-A15-18 *4	18	457
FO-A15-21 *4	21	533
FO-A15-24 *4	24	610




Remote control unit - 1 door

Model no.	Width	
	inches	mm
FO-A24-18	18	457
FO-A24-24	24	610

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.




1 door sink unit
- false panel above

FO-A15A-15 *4	15	381
FO-A15A-18 *4	18	457
FO-A15A-21 *4	21	533
FO-A15A-24 *4	24	610

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.




Remote control unit - 2 doors

FO-A24-30	30	762
FO-A24-36	36	914
FO-A24-48	48	1219


Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.



2 door sink unit

FO-A16-30 *4	30	762
FO-A16-36 *4	36	914
FO-A16-42 *4	42	1067
FO-A16-48 *4	48	1219




1 door bin unit

FO-A28-18	18	457
FO-A28-24	24	316

Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SSB" for stainless steel type 304;
 - For leak proof bin add the letters "LP".




2 door sink unit
- false panel above

FO-A16A-30 *4	30	762
FO-A16A-36 *4	36	914
FO-A16A-42 *4	42	1067
FO-A16A-48 *4	48	1219

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.




2 door bin unit

FO-A29-36	36	914
FO-A29-48	48	1219


Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SSB" for stainless steel type 304;
 - For leak proof bin add the letters "LP".




Island end sink unit

FO-A17-52 *4	52	1321
FO-A17-58 *4	58	1473



Open unit


FO-A35-15	15	381
FO-A35-18	18	457
FO-A35-21	21	533
FO-A35-24	24	610
FO-A35-30	30	762
FO-A35-36	36	914
FO-A35-42	42	1067
FO-A35-48	48	1219



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Sitting Height Base Cabinets




Model no.	Width	
	inches	mm
FO-B1-12-18	12	305
FO-B1-12-21	12	305
FO-B1-15	15	381
FO-B1-18	18	457
FO-B1-21	21	533
FO-B1-24	24	610

1 door unit



FO-B2-30	30	762
FO-B2-36	36	914
FO-B2-42	42	1067
FO-B2-48	48	1219

2 door unit




FO-B4-15-1	15	381
FO-B4-18-1	18	457
FO-B4-21-1	21	533
FO-B4-24-1	24	610

1 drawer & 1 door unit




FO-B4-15-2	15	381
FO-B4-18-2	18	457
FO-B4-21-2	21	533
FO-B4-24-2	24	610

2 half drawer & 1 door unit



FO-B5-30-1 *1	30	762
FO-B5-36-1 *1	36	914
FO-B5-42-1 *1*3	42	1067
FO-B5-48-1 *1*3	48	1219

1 drawer & 2 door unit




Model no.	Width	
	inches	mm
FO-B5-30-2	30	762
FO-B5-36-2	36	914
FO-B5-42-2	42	1067
FO-B5-48-2	48	1219

2 drawer & 2 door unit




FO-B5-30-4	30	762
FO-B5-36-4	36	914
FO-B5-42-4	42	1067
FO-B5-48-4	48	1219

4 half drawer & 2 door unit




FO-B6-30-3	30	762
FO-B6-36-3	36	914
FO-B6-42-3	42	1067
FO-B6-48-3	48	1219

3 drawer & 1 door unit




FO-B6-30-4	30	762
FO-B6-36-4	36	914
FO-B6-42-4	42	1067
FO-B6-48-4	48	1219

2 half drawer, 2 drawer & 1 door unit



FO-B6-30-5	30	762
FO-B6-36-5	36	914
FO-B6-42-5	42	1067
FO-B6-48-5	48	1219


4 half drawer, 1 drawer & 1 door unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.




Model no.	Width	
	inches	mm
FO-B6-30-6	30	762
FO-B6-36-6	36	914
FO-B6-42-6	42	1067
FO-B6-48-6	48	1219

6 half drawer & 1 door unit




FO-B6A-30-4	30	762
FO-B6A-36-4	36	914
FO-B6A-42-4	42	1067
FO-B6A-48-4	48	1219




Model no.	Width	
	inches	mm
FO-B8-30-4 *1	30	762
FO-B8-36-4 *1	36	914
FO-B8-42-4 *1*3	42	1067
FO-B8-48-4 *1*3	48	1219

4 drawer unit




FO-B8-30-6 *1	30	762
FO-B8-36-6 *1	36	914
FO-B8-42-6 *1*3	42	1067
FO-B8-48-6 *1*3	48	1219

4 drawer & 1 door unit




FO-B7-15-3	15	381
FO-B7-18-3	18	457
FO-B7-21-3	21	533
FO-B7-24-3	24	610
FO-B7-30-3 *1	30	762
FO-B7-36-3 *1	36	914
FO-B7-42-3 *1*3	42	1067
FO-B7-48-3 *1*3	48	1219

3 drawer unit



FO-B7-15-4	15	381
FO-B7-18-4	18	457
FO-B7-21-4	21	533
FO-B7-24-4	24	610
FO-B7-30-4 *1	30	762
FO-B7-36-4 *1	36	914
FO-B7-42-4 *1*3	42	1067
FO-B7-48-4 *1*3	48	1219


2 half drawer & 2 drawer unit



FO-B7-15-6	15	381
FO-B7-18-6	18	457
FO-B7-21-6	21	533
FO-B7-24-6	24	610
FO-B7-30-6 *1	30	762
FO-B7-36-6 *1	36	914
FO-B7-42-6 *1*3	42	1067
FO-B7-48-6 *1*3	48	1219


6 half drawer unit

4 half drawer & 2 drawer unit



FO-B10-30-6	30	762
FO-B10-36-6	36	914
FO-B10-42-6	42	1067
FO-B10-48-6	48	1219

6 drawer unit




FO-B10-30-8	30	762
FO-B10-36-8	36	914
FO-B10-42-8	42	1067
FO-B10-48-8	48	1219

4 half drawer & 4 drawer unit



FO-B10-30-12	30	762
FO-B10-36-12	36	914
FO-B10-42-12	42	1067
FO-B10-48-12	48	1219

12 half drawer unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Sitting Height Base Cabinets




Model no.	Width	
	inches	mm
FO-B11-15-2 *2	15	381
FO-B11-18-2 *2	18	457
FO-B11-21-2 *2	21	533
FO-B11-24-2 *2	24	610

1 drawer & 1 file drawer unit



FO-B11-15-3 *2	15	381
FO-B11-18-3 *2	18	457
FO-B11-21-3 *2	21	533
FO-B11-24-3 *2	24	610

2 half drawer & 1 file drawer unit




FO-B11-30-4 *2	30	762
FO-B11-36-4 *2	36	914
FO-B11-42-4 *2	42	1067
FO-B11-48-4 *2	48	1219

2 drawer & 2 file drawer unit




FO-B11-30-6 *2	30	762
FO-B11-36-6 *2	36	914
FO-B11-42-6 *2	42	1067
FO-B11-48-6 *2	48	1219

4 half drawer & 2 file drawer unit




FO-B20-30	30	762
FO-B20-36	36	914
FO-B20-42	42	1067
FO-B20-48	48	1219

Corner unit



Model no.	Width	
	inches	mm
FO-B21-30	30	762
FO-B21-36	36	914
FO-B21-42	42	1067
FO-B21-48	48	1219

Corner unit




FO-B28-18	18	457
FO-B28-24	24	316

Description


- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

1 door bin unit




FO-B35-15	15	381
FO-B35-18	18	457
FO-B35-21	21	533
FO-B35-24	24	610
FO-B35-30	30	762
FO-B35-36	36	914
FO-B35-42	42	1067
FO-B35-48	48	1219

Open unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FO-B15-15 *4	15	381
FO-B15-18 *4	18	457
FO-B15-21 *4	21	533
FO-B15-24 *4	24	610

1 door sink unit




FO-B16-30 *4	30	762
FO-B16-36 *4	36	914
FO-B16-42 *4	42	1067
FO-B16-48 *4	48	1219

2 door sink unit



FO-B17-52 *4	52	1321
FO-B17-58 *4	58	1473


Island end sink unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

ADA Height Base Cabinets




Model no.	Width	
	inches	mm
FO-ADA1-12-18	12	305
FO-ADA1-12-21	12	305
FO-ADA1-15	15	381
FO-ADA1-18	18	457
FO-ADA1-21	21	533
FO-ADA1-24	24	610

1 door unit




FO-ADA2-30	30	762
FO-ADA2-36	36	914
FO-ADA2-42	42	1067
FO-ADA2-48	48	1219

2 door unit




FO-ADA4-15-1	15	381
FO-ADA4-18-1	18	457
FO-ADA4-21-1	21	533
FO-ADA4-24-1	24	610

1 top drawer & 1 door unit




FO-ADA5-30-1	30	762
FO-ADA5-36-1	36	914
FO-ADA5-42-1 *3	42	1067
FO-ADA5-48-1 *3	48	1219

1 top drawer & 2 door unit




FO-ADA5-30-2	30	762
FO-ADA5-36-2	36	914
FO-ADA5-42-2	42	1067
FO-ADA5-48-2	48	1219

2 top drawer & 2 door unit




Model no.	Width	
	inches	mm
FO-ADA6-30-4	30	762
FO-ADA6-36-4	36	914
FO-ADA6-42-4	42	1067
FO-ADA6-48-4	48	1219

3 drawer, 1 top drawer & 1 door unit




FO-ADA6A-30-5	30	762
FO-ADA6A-36-5	36	914
FO-ADA6A-42-5	42	1067
FO-ADA6A-48-5	48	1219

3 drawer, 2 top drawer & 1 door unit




FO-ADA6-30-7	30	762
FO-ADA6-36-7	36	914
FO-ADA6-42-7	42	1067
FO-ADA6-48-7	48	1219

6 half drawer, 1 top drawer & 1 door unit




FO-ADA7-15-4	15	381
FO-ADA7-18-4	18	457
FO-ADA7-21-4	21	533
FO-ADA7-24-4	24	610
FO-ADA7-30-4 *1	30	762
FO-ADA7-36-4 *1	36	914
FO-ADA7-42-4 *1*3	42	1067
FO-ADA7-48-4 *1*3	48	1219

3 drawer & 1 top drawer unit



FO-ADA8-30-5	30	762
FO-ADA8-36-5	36	914
FO-ADA8-42-5 *3	42	1067
FO-ADA8-48-5 *3	48	1219


3 drawer & 2 top drawer unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.




Model no.	Width	
	inches	mm
FO-ADA9-30-6	30	762
FO-ADA9-36-6	36	914
FO-ADA9-42-6 *3	42	1067
FO-ADA9-48-6 *3	48	1219

4 drawer & 2 top drawer unit




FO-ADA10-30-8	30	762
FO-ADA10-36-8	36	914
FO-ADA10-42-8	42	1067
FO-ADA10-48-8	48	1219

6 drawer & 2 top drawer unit




FO-ADA11-15-3 *2	15	381
FO-ADA11-18-3 *2	18	457
FO-ADA11-21-3 *2	21	533
FO-ADA11-24-3 *2	24	610

1 drawer, 1 top drawer & 1 file drawer unit




FO-ADA11-30-6 *2	30	762
FO-ADA11-36-6 *2	36	914
FO-ADA11-42-6 *2	42	1067
FO-ADA11-48-6 *2	48	1219

2 drawer, 2 top drawer & 2 file drawer unit



FO-ADA12-30-3 *2	30	762
FO-ADA12-36-3 *2	36	914
FO-ADA12-42-3 *2	42	1067
FO-ADA12-48-3 *2	48	1219

1 drawer, 1 top drawer, 1 file drawer & 1 door unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

ADA Height Base Cabinets




Model no.	Width	
	inches	mm
FO-ADA15-15 *4	15	381
FO-ADA15-18 *4	18	457
FO-ADA15-21 *4	21	533
FO-ADA15-24 *4	24	610

1 door sink unit




FO-ADA15A-15 *4	15	381
FO-ADA15A-18 *4	18	457
FO-ADA15A-21 *4	21	533
FO-ADA15A-24 *4	24	610

1 door sink unit - false panel above




FO-ADA16-30 *4	30	762
FO-ADA16-36 *4	36	914
FO-ADA16-42 *4	42	1067
FO-ADA16-48 *4	48	1219

2 door sink unit




FO-ADA16A-30 *4	30	762
FO-ADA16A-36 *4	36	914
FO-ADA16A-42 *4	42	1067
FO-ADA16A-48 *4	48	1219

2 door sink unit - false panel above




FO-ADA17-52 *4	52	1321
FO-ADA17-58 *4	58	1473

Island end sink unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FO-ADA20-30	30	762
FO-ADA20-36	36	914
FO-ADA20-42	42	1067
FO-ADA20-48	48	1219

Corner unit




FO-ADA21-30	30	762
FO-ADA21-36	36	914
FO-ADA21-42	42	1067
FO-ADA21-48	48	1219

Corner unit



FO-ADA24-18	18	457
FO-ADA24-24	24	610

Remote control unit - 1 door




FO-ADA24-30	30	762
FO-ADA24-36	36	914
FO-ADA24-48	48	1219

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

Remote control unit - 2 doors




FO-ADA28-18	18	457
FO-ADA28-24	24	316

Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

1 door bin unit




Model no.	Width	
	inches	mm
FO-ADA29-36	36	914
FO-ADA29-48	48	1219

Description


- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

2 door bin unit



FO-ADA35-15	15	381
FO-ADA35-18	18	457
FO-ADA35-21	21	533
FO-ADA35-24	24	610
FO-ADA35-30	30	762
FO-ADA35-36	36	914
FO-ADA35-42	42	1067
FO-ADA35-48	48	1219


Open unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

24-3/4" High Suspended Cabinets




Model no.	Width	
	inches	mm
FO-OB1-15	15	381
FO-OB1-18	18	914
FO-OB1-21	21	1067
FO-OB1-24	24	1219




Model no.	Width	
	inches	mm
FO-OB10-30-6	30	752
FO-OB10-36-6	36	914
FO-OB10-42-6	42	1067
FO-OB10-48-6	48	1219

Suspended 1 Hinged Door Unit




FO-OB2-30	30	752
FO-OB2-36	36	914
FO-OB2-42	42	1067
FO-OB2-48	48	1219

Suspended 6 Drawer Unit




FO-OB10-30-8	30	752
FO-OB10-36-8	36	914
FO-OB10-42-8	42	1067
FO-OB10-48-8	48	1219

Suspended 2 Hinged Door unit




FO-OB4-15-1 *2	15	381
FO-OB4-18-1 *2	18	457
FO-OB4-21-1 *2	21	533
FO-OB4-24-1 *2	24	610

Suspended 4 Half Drawer & 4 Regular Drawer unit




FO-OB11-15-2 *2	15	381
FO-OB11-18-2 *2	18	457
FO-OB11-21-2 *2	21	533
FO-OB11-24-2 *2	24	610

Suspended 1 Drawer & 1 Hinged Door unit




FO-OB5-30-2	30	752
FO-OB5-36-2	36	914
FO-OB5-42-2	42	1067
FO-OB5-48-2	48	1219

Suspended 1 Drawer & 1 File Drawer unit




FO-OB11-15-3 *2	15	381
FO-OB11-18-3 *2	18	457
FO-OB11-21-3 *2	21	533
FO-OB11-24-3 *2	24	610

Suspended 2 Drawer & 2 Hinged Door unit



FO-OB7-15-3	15	381
FO-OB7-18-3	18	457
FO-OB7-21-3	21	533
FO-OB7-24-3	24	610
FO-OB7-30-3	30	752
FO-OB7-36-3	36	914
FO-OB7-42-3	42	1067
FO-OB7-48-3	48	1219

Suspended 2 Half Drawer & 1 File Drawer unit




FO-OB35-15	15	381
FO-OB35-18	18	457
FO-OB35-21	21	533
FO-OB35-24	24	610
FO-OB35-30	30	752
FO-OB35-36	36	914
FO-OB35-42	42	1067
FO-OB35-48	48	1219

Suspended 3 Drawer unit




FO-OB7-15-4	15	381
FO-OB7-18-4	18	457
FO-OB7-21-4	21	533
FO-OB7-24-4	24	610
FO-OB7-30-4	30	752
FO-OB7-36-4	36	914
FO-OB7-42-4	42	1067
FO-OB7-48-4	48	1214

Suspended 2 Half Drawer & 2 Drawer unit



Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FO-OC1-15	15	381
FO-OC1-18	18	457
FO-OC1-21	21	533
FO-OC1-24	24	610




Model no.	Width	
	inches	mm
FO-OC10-30-6	30	762
FO-OC10-36-6	36	914
FO-OC10-42-6	42	1067
FO-OC10-48-6	48	1219

Suspended 1 Hinged Door unit




FO-OC2-30	30	762
FO-OC2-36	36	914
FO-OC2-42	42	1067
FO-OC2-48	48	1219

Suspended 2 Half Drawer & 4 Drawer unit




FO-OC11-15-2 *2	15	381
FO-OC11-18-2 *2	18	457
FO-OC11-21-2 *2	21	533
FO-OC11-24-2 *2	24	610

Suspended 2 Hinged Door unit




FO-OC7-15-3	15	381
FO-OC7-18-3	18	457
FO-OC7-21-3	21	533
FO-OC7-24-3	24	610

Suspended 1 Half Drawer + 1 File Drawer unit




FO-OC11-30-4 *2	30	762
FO-OC11-36-4 *2	36	914
FO-OC11-42-4 *2	42	1067
FO-OC11-48-4 *2	48	1219

Suspended 1 Half Drawer & 2 Drawer unit



FO-OC10-15-4	15	381
FO-OC10-18-4	18	457
FO-OC10-21-4	21	533
FO-OC10-24-4	24	610


Suspended 2 Half Drawer + 2 File Drawer unit



FO-OC35-15	15	381
FO-OC35-18	18	457
FO-OC35-21	21	533
FO-OC35-24	24	610
FO-OC35-30	30	762
FO-OC35-36	36	914
FO-OC35-42	42	1067
FO-OC35-48	48	1219

Suspended 3 Half Drawer & 1 Drawer unit

Suspended Open unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

24" High Wall cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FO-E10-24-15-14	15	381
FO-E10-24-18-14	18	457
FO-E10-24-21-14	21	533
FO-E10-24-24-14	24	610
FO-E10-24-30-14	30	762
FO-E10-24-36-14	36	914
FO-E10-24-42-14	42	1067
FO-E10-24-48-14	48	1219
FO-E10-24-15-18	15	381
FO-E10-24-18-18	18	457
FO-E10-24-21-18	21	533
FO-E10-24-24-18	24	610
FO-E10-24-30-18	30	762
FO-E10-24-36-18	36	914
FO-E10-24-42-18	42	1067
FO-E10-24-48-18	48	1219




2 hinged door unit

Model no.	Width	
	inches	mm
FO-E20-24-30-14	30	762
FO-E20-24-36-14	36	914
FO-E20-24-42-14	42	1067
FO-E20-24-48-14	48	1219
FO-E20-24-30-18	30	762
FO-E20-24-36-18	36	914
FO-E20-24-42-18	42	1067
FO-E20-24-48-18	48	1219




2 glazed hinged door unit

FO-E30-24-30-14	30	762
FO-E30-24-36-14	36	914
FO-E30-24-42-14	42	1067
FO-E30-24-48-14	48	1219
FO-E30-24-30-18	30	762
FO-E30-24-36-18	36	914
FO-E30-24-42-18	42	1067
FO-E30-24-48-18	48	1219




1 hinged door unit

FO-E20-24-15-14	15	381
FO-E20-24-18-14	18	457
FO-E20-24-21-14	21	533
FO-E20-24-24-14	24	610
FO-E20-24-15-18	15	381
FO-E20-24-18-18	18	457
FO-E20-24-21-18	21	533
FO-E20-24-24-18	24	610




1 glazed hinged door unit

FO-E30-24-15-14	15	381
FO-E30-24-18-14	18	457
FO-E30-24-21-14	21	533
FO-E30-24-24-14	24	610
FO-E30-24-15-18	15	381
FO-E30-24-18-18	18	457
FO-E30-24-21-18	21	533
FO-E30-24-24-18	24	610




Corner unit with door

FO-E22-24-30-14	30	762
FO-E22-24-36-14	36	914
FO-E22-24-42-14	42	1067
FO-E22-24-48-14	48	1219
FO-E22-24-30-18	30	762
FO-E22-24-36-18	36	914
FO-E22-24-42-18	42	1067
FO-E22-24-48-18	48	1219




Corner unit with door

FO-E32-24-30-14	30	762
FO-E32-24-36-14	36	914
FO-E32-24-42-14	42	1067
FO-E32-24-48-14	48	1219
FO-E32-24-30-18	30	762
FO-E32-24-36-18	36	914
FO-E32-24-42-18	42	1067
FO-E32-24-48-18	48	1219




14" depth units




18" depth units

Reverse models are also available by adding "RV" to model number.



Open shelf unit

Model no.	Width	
	inches	mm
FO-E10-30-15-14	15	381
FO-E10-30-18-14	18	457
FO-E10-30-21-14	21	533
FO-E10-30-24-14	24	610
FO-E10-30-30-14	30	762
FO-E10-30-36-14	36	914
FO-E10-30-42-14	42	1067
FO-E10-30-48-14	48	1219
FO-E10-30-15-18	15	381
FO-E10-30-18-18	18	457
FO-E10-30-21-18	21	533
FO-E10-30-24-18	24	610
FO-E10-30-30-18	30	762
FO-E10-30-36-18	36	914
FO-E10-30-42-18	42	1067
FO-E10-30-48-18	48	1219




2 glazed hinged door unit

Model no.	Width	
	inches	mm
FO-E30-30-30-14	30	762
FO-E30-30-36-14	36	914
FO-E30-30-42-14	42	1067
FO-E30-30-48-14	48	1219
FO-E30-30-30-18	30	762
FO-E30-30-36-18	36	914
FO-E30-30-42-18	42	1067
FO-E30-30-48-18	48	1219


Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




1 hinged door unit

FO-E20-30-15-14	15	381
FO-E20-30-18-14	18	457
FO-E20-30-21-14	21	533
FO-E20-30-24-14	24	610
FO-E20-30-15-18	15	381
FO-E20-30-18-18	18	457
FO-E20-30-21-18	21	533
FO-E20-30-24-18	24	610




1 glazed hinged door unit

FO-E30-30-15-14	15	381
FO-E30-30-18-14	18	457
FO-E30-30-21-14	21	533
FO-E30-30-24-14	24	610
FO-E30-30-15-18	15	381
FO-E30-30-18-18	18	457
FO-E30-30-21-18	21	533
FO-E30-30-24-18	24	610




Corner unit with door

FO-E22-30-30-14	30	762
FO-E22-30-36-14	36	914
FO-E22-30-42-14	42	1067
FO-E22-30-48-14	48	1219
FO-E22-30-30-18	30	762
FO-E22-30-36-18	36	914
FO-E22-30-42-18	42	1067
FO-E22-30-48-18	48	1219




2 hinged door unit

FO-E20-30-30-14	30	762
FO-E20-30-36-14	36	914
FO-E20-30-42-14	42	1067
FO-E20-30-48-14	48	1219
FO-E20-30-30-18	30	762
FO-E20-30-36-18	36	914
FO-E20-30-42-18	42	1067
FO-E20-30-48-18	48	1219




Corner unit with glazed door

FO-E32-30-30-14	30	762
FO-E32-30-36-14	36	914
FO-E32-30-42-14	42	1067
FO-E32-30-48-14	48	1219
FO-E32-30-30-18	30	762
FO-E32-30-36-18	36	914
FO-E32-30-42-18	42	1067
FO-E32-30-48-18	48	1219




14" depth units



18" depth units


Reverse models are also available by adding "RV" to model number.

36" High Wall Cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FO-E10-36-15-14	15	381
FO-E10-36-18-14	18	457
FO-E10-36-21-14	21	533
FO-E10-36-24-14	24	610
FO-E10-36-30-14	30	762
FO-E10-36-36-14	36	914
FO-E10-36-42-14	42	1067
FO-E10-36-48-18	48	1219
FO-E10-36-15-18	15	381
FO-E10-36-18-18	18	457
FO-E10-36-21-18	21	533
FO-E10-36-24-18	24	610
FO-E10-36-30-18	30	762
FO-E10-36-36-18	36	914
FO-E10-36-42-18	42	1067
FO-E10-36-48-18	48	1219



2 solid hinged door unit

Model no.	Width	
	inches	mm
FO-E20-36-30-14	30	762
FO-E20-36-36-14	36	914
FO-E20-36-42-14	42	1067
FO-E20-36-48-14	48	1219
FO-E20-36-30-18	30	762
FO-E20-36-36-18	36	914
FO-E20-36-42-18	42	1067
FO-E20-36-48-18	48	1219




2 glazed hinged door unit

FO-E30-36-30-14	30	762
FO-E30-36-36-14	36	914
FO-E30-36-42-14	42	1067
FO-E30-36-48-14	48	1219
FO-E30-36-30-18	30	762
FO-E30-36-36-18	36	914
FO-E30-36-42-18	42	1067
FO-E30-36-48-18	48	1219


Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




1 hinged door unit

FO-E20-36-15-14	15	381
FO-E20-36-18-14	18	457
FO-E20-36-21-14	21	533
FO-E20-36-24-14	24	610
FO-E20-36-15-18	15	381
FO-E20-36-18-18	18	457
FO-E20-36-21-18	21	533
FO-E20-36-24-18	24	610




1 glazed hinged door unit

FO-E30-36-15-14	15	381
FO-E30-36-18-14	18	457
FO-E30-36-21-14	21	533
FO-E30-36-24-14	24	610
FO-E30-36-15-18	15	381
FO-E30-36-18-18	18	457
FO-E30-36-21-18	21	533
FO-E30-36-24-18	24	610




Corner unit with door

FO-E22-36-30-14	30	762
FO-E22-36-36-14	36	914
FO-E22-36-42-14	42	1067
FO-E22-36-48-14	48	1219
FO-E22-36-30-18	30	762
FO-E22-36-36-18	36	914
FO-E22-36-42-18	42	1067
FO-E22-36-48-18	48	1219




Corner unit with door


FO-E32-36-30-14	30	762
FO-E32-36-36-14	36	914
FO-E32-36-42-14	42	1067
FO-E32-36-48-14	48	1219
FO-E32-36-30-18	30	762
FO-E32-36-36-18	36	914
FO-E32-36-42-18	42	1067
FO-E32-36-48-18	48	1219



14" depth units




18" depth units




Open shelf unit

Model no.	Width	
	inches	mm
FO-E10-48-15-14	15	381
FO-E10-48-18-14	18	457
FO-E10-48-21-14	21	533
FO-E10-48-24-14	24	610
FO-E10-48-30-14	30	762
FO-E10-48-36-14	36	914
FO-E10-48-42-14	42	1067
FO-E10-48-48-14	48	1219
FO-E10-48-15-18	15	381
FO-E10-48-18-18	18	457
FO-E10-48-21-18	21	533
FO-E10-48-24-18	24	610
FO-E10-48-30-18	30	762
FO-E10-48-36-18	36	914
FO-E10-48-42-18	42	1067
FO-E10-48-48-18	48	1219




2 hinged door unit

Model no.	Width	
	inches	mm
FO-E20-48-30-14	30	762
FO-E20-48-36-14	36	914
FO-E20-48-42-14	42	1067
FO-E20-48-48-14	48	1219
FO-E20-48-30-18	30	762
FO-E20-48-36-18	36	914
FO-E20-48-42-18	42	1067
FO-E20-48-48-18	48	1219




2 glazed hinged door unit

FO-E30-48-30-14	30	762
FO-E30-48-36-14	36	914
FO-E30-48-42-14	42	1067
FO-E30-48-48-14	48	1219
FO-E30-48-30-18	30	762
FO-E30-48-36-18	36	914
FO-E30-48-42-18	42	1067
FO-E30-48-48-18	48	1219



1 hinged door unit

FO-E20-48-15-14	15	381
FO-E20-48-18-14	18	457
FO-E20-48-21-14	21	533
FO-E20-48-24-14	24	610
FO-E20-48-15-18	15	381
FO-E20-48-18-18	18	457
FO-E20-48-21-18	21	533
FO-E20-48-24-18	24	610




1 glazed hinged door unit

FO-E30-48-15-14	15	381
FO-E30-48-18-14	18	457
FO-E30-48-21-14	21	533
FO-E30-48-24-14	24	610
FO-E30-48-15-18	15	381
FO-E30-48-18-18	18	457
FO-E30-48-21-18	21	533
FO-E30-48-24-18	24	610


Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




Corner unit - with door 48"

FO-E22-48-30-14	30	762
FO-E22-48-36-14	36	914
FO-E22-48-42-14	42	1067
FO-E22-48-48-14	48	1219
FO-E22-48-30-18	30	762
FO-E22-48-36-18	36	914
FO-E22-48-42-18	42	1067
FO-E22-48-48-18	48	1219



Corner unit - with glazed door 48"

FO-E32-48-30-14	30	762
FO-E32-48-36-14	36	914
FO-E32-48-42-14	42	1067
FO-E32-48-48-14	48	1219
FO-E32-48-30-18	30	762
FO-E32-48-36-18	36	914
FO-E32-48-42-18	42	1067
FO-E32-48-48-18	48	1219




14" depth units



18" depth units


Reverse models are also available by adding "RV" to model number.

84" High Tall Cabinets




Open shelf tall cab

Model no.	Width	
	inches	mm
FO-F10-84-15-14	15	381
FO-F10-84-18-14	18	457
FO-F10-84-21-14	21	533
FO-F10-84-24-14	24	610
FO-F10-84-30-14	30	762
FO-F10-84-36-14	36	914
FO-F10-84-42-14	42	1067
FO-F10-84-48-14	48	1219
FO-F10-84-15-18	15	381
FO-F10-84-18-18	18	457
FO-F10-84-21-18	21	533
FO-F10-84-24-18	24	610
FO-F10-84-30-18	30	762
FO-F10-84-36-18	36	914
FO-F10-84-42-18	42	1067
FO-F10-84-48-18	48	1219
FO-F10-84-15-22	15	381
FO-F10-84-18-22	18	457
FO-F10-84-21-22	21	533
FO-F10-84-24-22	24	610
FO-F10-84-30-22	30	762
FO-F10-84-36-22	36	914
FO-F10-84-42-22	42	1067
FO-F10-84-48-22	48	1219




Model no.	Width	
	inches	mm
FO-F20-84-30-14	30	762
FO-F20-84-36-14	36	914
FO-F20-84-42-14	42	1067
FO-F20-84-48-14	48	1219
FO-F20-84-30-18	30	762
FO-F20-84-36-18	36	914
FO-F20-84-42-18	42	1067
FO-F20-84-48-18	48	1219
FO-F20-84-30-22	30	762
FO-F20-84-36-22	36	914
FO-F20-84-42-22	42	1067
FO-F20-84-48-22	48	1219

2 hinged door tall cab




FO-F21-84-15-14	15	381
FO-F21-84-18-14	18	457
FO-F21-84-21-14	21	533
FO-F21-84-24-14	24	610
FO-F21-84-15-18	15	381
FO-F21-84-18-18	18	457
FO-F21-84-21-18	21	533
FO-F21-84-24-18	24	610
FO-F21-84-15-22	15	381
FO-F21-84-18-22	18	457
FO-F21-84-21-22	21	533
FO-F21-84-24-22	24	610

2 hinged door tall cab



FO-F20-84-15-14	15	381
FO-F20-84-18-14	18	457
FO-F20-84-21-14	21	533
FO-F20-84-24-14	24	610
FO-F20-84-15-18	15	381
FO-F20-84-18-18	18	457
FO-F20-84-21-18	21	533
FO-F20-84-24-18	24	610
FO-F20-84-15-22	15	381
FO-F20-84-18-22	18	457
FO-F20-84-21-22	21	533
FO-F20-84-24-22	24	610

1 hinged door tall cab




FO-F21-84-30-14	30	762
FO-F21-84-36-14	36	914
FO-F21-84-42-14	42	1067
FO-F21-84-48-14	48	1219
FO-F21-84-30-18	30	762
FO-F21-84-36-18	36	914
FO-F21-84-42-18	42	1067
FO-F21-84-48-18	48	1219
FO-F21-84-30-22	30	762
FO-F21-84-36-22	36	914
FO-F21-84-42-22	42	1067
FO-F21-84-48-22	48	1219

4 hinged door tall cab




Tall cabinet available depths of 14", 18" and 22"




Model no.	Width	
	inches	mm
FO-F30-84-15-14	15	381
FO-F30-84-18-14	18	457
FO-F30-84-21-14	21	533
FO-F30-84-24-14	24	610
FO-F30-84-15-18	15	381
FO-F30-84-18-18	18	457
FO-F30-84-21-18	21	533
FO-F30-84-24-18	24	610
FO-F30-84-15-22	15	381
FO-F30-84-18-22	18	457
FO-F30-84-21-22	21	533
FO-F30-84-24-22	24	610

1 glazed hinged door tall cab




Model no.	Width	
	inches	mm
FO-F31-84-15-14	15	381
FO-F31-84-18-14	18	457
FO-F31-84-21-14	21	533
FO-F31-84-24-14	24	610
FO-F31-84-15-18	15	381
FO-F31-84-18-18	18	457
FO-F31-84-21-18	21	533
FO-F31-84-24-18	24	610
FO-F31-84-15-22	15	381
FO-F31-84-18-22	18	457
FO-F31-84-21-22	21	533
FO-F31-84-24-22	24	610

2 glazed hinged door tall cab



FO-F30-84-30-14	30	762
FO-F30-84-36-14	36	914
FO-F30-84-42-14	42	1067
FO-F30-84-48-14	48	1219
FO-F30-84-30-18	30	762
FO-F30-84-36-18	36	914
FO-F30-84-42-18	42	1067
FO-F30-84-48-18	48	1219
FO-F30-84-30-22	30	762
FO-F30-84-36-22	36	914
FO-F30-84-42-22	42	1067
FO-F30-84-48-22	48	1219

2 glazed hinged door tall cab



FO-F31-84-30-14	30	762
FO-F31-84-36-14	36	914
FO-F31-84-42-14	42	1067
FO-F31-84-48-14	48	1219
FO-F31-84-30-18	30	762
FO-F31-84-36-18	36	914
FO-F31-84-42-18	42	1067
FO-F31-84-48-18	48	1219
FO-F31-84-30-22	30	762
FO-F31-84-36-22	36	914
FO-F31-84-42-22	42	1067
FO-F31-84-48-22	48	1219

4 glazed hinged door tall cab

Description


- Framed glass doors to be supplied using 3 mm float glass



Tall cabinet: available depths of 14", 18" and 22"

1 CABINET PULLS**Standard — Brushed nickel wire pull 4" c/c**


- WP ☐ Stainless steel wire pull 4" c/c
 PP ☐ Black recessed 4" x 2" pull
 AP ☐ Aluminum recessed 4" x 2" pull

**2 DRAWER SLIDE ARMS****Standard — 75% opening self closing slide arms on all regular drawers – SEFA 8M – 150 lbs****Standard — 100% full extension opening on file drawers – SEFA 8M – 200 lbs****Standard — 100% full extension opening on all drawers larger than 36" – SEFA 8M – 200 lbs**

- TZ ☐ 100% full extension opening on all project drawers – SEFA 8M – 100 lbs
 TY ☐ 100% full extension opening with self closing option on all project drawers – SEFA 8M – 100 lbs
 TX ☐ 100% full extension opening on all project drawers – SEFA 8M – 200 lbs

3 CABINET HINGES**Standard — knuckles nickel plated finish**

- FH ☐ 5-knuckles stainless steel finish
 MH ☐ Polyamide mounted nickel plated steel MB-20

**4 LOCKS****Standard — None**

- KA ☐ 8-disc chrome plated with removable core on all doors and all drawers
 KT ☐ 8-disc chrome plated with removable core on all drawers
 KP ☐ 8-disc chrome plated with removable core on all doors and drawers
 K1 ☐ 8-disc chrome plated with removable core on top drawers only
 KF ☐ 8-disc chrome plated with removable core on file drawers only
 KC ☐ Chrome plated padlock bracket on all doors and all drawers
 XI ☐ Interlocking system that permits only one drawer opening at the time (with full extension arms only) – standard on mobile cabinets
 TP ☐ Theft security panels between drawers
 L1 ☐ All locks different with a project master key
 L2 ☐ All locks the same per cabinet with a project master key
 L3 ☐ All locks the same per room with a room master key
 L4 ☐ All locks the same

5 DOOR CATCHES**Standard — Spring loaded roller door catch****6 CARD HOLDERS****Standard — None**

- CH ☐ Aluminum 4" x 2" card holder on all drawers
 CD ☐ Aluminum 4" x 2" card holder on all doors



7 CABINET SHELF

Standard — One (1) full depth shelf, adjustable on ½" increment

- T1 ☐ Additional shelf per cabinet
- P1 ☐ Pull out shelf for cabinets 15", 18", 21" or 24" only
- PR ☐ Pull out shelf on right cabinet side
- PL ☐ Pull out shelf on left cabinet side
- S2 ☐ Split shelves (2 pieces), adjustable on ½" increment
- FS ☐ Cabinet full extension shelf supplied with full extension slide arms – 200 lbs

8 FINISHED PANELS

Standard — Removable access panel on cupboard cabinets only (open back on drawer cab)

- BP ☐ Removable access panel on drawer units
- FO ☐ Painted steel flat finished back welded to cabinet
- FT ☐ Painted steel flat finished top and back panels welded to cabinet

9 CABINET BASE OPTIONS

Standard — 3" x 4" high integral front toe kick with four 3/8"x 2 ½" leveler glides

- PN ☐ Cabinet levelers supplied with grey rubber cap
- SS ☐ Four swivel casters 4" height, rubber, with brakes on two front casters – 275 lbs
- WB ☐ Cabinet supplied with no toe kick (4" shorter)
- RT ☐ Cabinet supplied with a right toe kick 3" x 4" high – not available on tall cabinets
- LT ☐ Cabinet supplied with a left toe kick 3" x 4" high – not available on tall cabinets
- VT ☐ Cabinet supplied with a venting grill on the front toe kick

10 VARIOUS OTHER CABINET OPTIONS

- RV ☐ Reverse door/drawer arrangement
- SK ☐ Sink support consisting of 2 horizontal struts with 4 rod hangers
- LG ☐ 6mm Laminated glass in lieu of standard 3mm float glass
- TG ☐ 6mm Tempered glass in lieu of standard 3mm float glass
- SB ☐ Stainless steel bin insert into bin cabinets
- GB ☐ Galvanized bin insert into bin cabinets
- LP ☐ Leak proof bin insert into bin cabinets
- WC ☐ Grill(s) incorporated into cabinet door(s)
- SS4 ☐ Same cabinet but supplied in stainless steel T.304
- SS6 ☐ Same cabinet but supplied in stainless steel T.316



BEDCOLAB LTD

**2305, Francis-Hugues Avenue
Laval, Quebec, Canada H7S 1N5**

**Phone: 514.384.2820 / 1.800.461.6414
Fax: 514.384.9210**

**E-mail: information@bedcolab.com
Web Site: www.bedcolab.com**



Authorized dealer:

--