





Features


Benefits

1	Superior front horizontal rail.	Additional rigidity. Impurity prevention.
2	Removable front center post.	Full flexibility for all door and drawer arrangements.
3	The front posts can accommodate varying arrangements of doors & drawers.	Flexibility even after installation.
4	Wrap-around cabinet shell construction.	Provides maximum strength and eliminates need for fixed finished back.
5	Removable back panel.	Allows accessibility to mechanical services.
6	Slotted rear post.	Support for slide arms and shelves. Full flexibility for future needs.
7	Cabinet body is fully enclosed.	Facilitating maintenance and cleaning.
8	Extended drawer depth.	Maximized storage capacity.
9	Insulated facades.	Operational noise reduced to a minimum.
10	Self-closing ball bearing slide arms.	Avoiding accidental spills
11	Threaded 5/16" leveler with large foot.	Solid, sturdy and prevents scratching.
12	Various standard handles.	Adhere to design intent.

Forte Inset Door/Drawer Combination Possibilities



Our cabinet design is engineered for interchangeability of the same drawer and door heights for all different cabinet types including the following drawer and door heights:



* The top drawer and full height door on an ADA cabinet are of different heights than for all other cabinets and will only fit at this position.

Reading the product identification numbers code

The Bedcolab product identification number is simple, making it easy to verify that you have ordered the correct cabinet. Each number or letter in the designation represents an aspect of the cabinet. Below are examples of the product identification numbers you will see in this catalogue along with an explanation of each designation:

Product Identification Code





FI	A	6A	30	5
Forte Inset Identification		Cabinet Height		
Cabinet Style		Quantity of Drawers		
Cabinet Style		Cabinet Width		
• A – standing • B – sitting • ADA – height		1 – 16		
• Drawer cabinet • Door cabinet • Combination door/drawer • Mixed combination door/drawer		15 in. 30 in. 18 in. 36 in. 21 in. 42 in. 24 in. 48 in.		

Base cabinets

FI	75	42	18
Forte Inset Identification		Cabinet Depth	
Cabinet Style		Cabinet Width	
• Single door • Double door • Glass door • Glazed door		14, 18, 22* (Tall cabinets)	
		15 in. 30 in. 18 in. 36 in. 21 in. 42 in. 24 in. 48 in.	


Wall and tall cabinets

Please refer to pages 38 and 39 of this Forte Inset catalogue for hardware options and optional cabinet features.




Typical Side Views

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FI-A1-12-18	12	305
FI-A1-12-21	12	305
FI-A1-15	15	381
FI-A1-18	18	457
FI-A1-21	21	533
FI-A1-24	24	610

1 door unit




Model no.	Width	
	inches	mm
FI-A5-30-1	30	762
FI-A5-36-1	36	914
FI-A5-42-1	42	1067
FI-A5-48-1	48	1219

1 drawer & 2 door unit




FI-A2-30	30	762
FI-A2-36	36	914
FI-A2-42	42	1067
FI-A2-48	48	1219

2 door unit




FI-A5-30-2	30	762
FI-A5-36-2	36	914
FI-A5-42-2	42	1067
FI-A5-48-2	48	1219

2 drawer & 2 door unit




FI-A3-30	30	762
FI-A3-36	36	914
FI-A3-42	42	1067
FI-A3-48	48	1219

2 sliding door unit




FI-A5-30-3	30	762
FI-A5-36-3	36	914
FI-A5-42-3	42	1067
FI-A5-48-3	48	1219

2 half drawer, 1 drawer & 2 door unit




FI-A4-15-1	15	381
FI-A4-18-1	18	457
FI-A4-21-1	21	533
FI-A4-24-1	24	610

1 drawer & 1 door unit




FI-A5-30-4	30	762
FI-A5-36-4	36	914
FI-A5-42-4	42	1067
FI-A5-48-4	48	1219

4 half drawer & 2 door unit




FI-A4-15-2	15	381
FI-A4-18-2	18	457
FI-A4-21-2	21	533
FI-A4-24-2	24	610

2 half drawer & 1 door unit



FI-A6-30-4	30	762
FI-A6-36-4	36	914
FI-A6-42-4	42	1067
FI-A6-48-4	48	1219


4 drawer & 1 door unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.




Model no.	Width	
	inches	mm
FI-A6-30-5	30	762
FI-A6-36-5	36	914
FI-A6-42-5	42	1067
FI-A6-48-5	48	1219

2 half drawer, 3 drawer & 1 door unit




Model no.	Width	
	inches	mm
FI-A7-15-4	15	381
FI-A7-18-4	18	457
FI-A7-21-4	21	533
FI-A7-24-4	24	610
FI-A7-30-4 *1	30	762
FI-A7-36-4 *1	36	914
FI-A7-42-4 *1*3	42	1067
FI-A7-48-4 *1*3	48	1219

4 drawer unit




FI-A6-30-7	30	762
FI-A6-36-7	36	914
FI-A6-42-7	42	1067
FI-A6-48-7	48	1219

6 half drawer, 1 drawer & 1 door unit




FI-A7-15-5	15	381
FI-A7-18-5	18	457
FI-A7-21-5	21	533
FI-A7-24-5	24	610
FI-A7-30-5 *1	30	762
FI-A7-36-5 *1	36	914
FI-A7-42-5 *1*3	42	1067
FI-A7-48-5 *1*3	48	1219

2 half drawer & 3 drawer unit




FI-A6A-30-5	30	762
FI-A6A-36-5	36	914
FI-A6A-42-5	42	1067
FI-A6A-48-5	48	1219

5 drawer & 1 door unit




FI-A7-15-6	15	381
FI-A7-18-6	18	457
FI-A7-21-6	21	533
FI-A7-24-6	24	610
FI-A7-30-6 *1	30	762
FI-A7-36-6 *1	36	914
FI-A7-42-6 *1*3	42	1067
FI-A7-48-6 *1*3	48	1219

4 half drawer & 2 drawer unit




FI-A6A-30-6	30	762
FI-A6A-36-6	36	914
FI-A6A-42-6	42	1067
FI-A6A-48-6	48	1219

2 half drawer, 4 drawer & 1 door unit




FI-A7-15-7	15	381
FI-A7-18-7	18	457
FI-A7-21-7	21	533
FI-A7-24-7	24	610
FI-A7-30-7 *1	30	762
FI-A7-36-7 *1	36	914
FI-A7-42-7 *1*3	42	1067
FI-A7-48-7 *1*3	48	1219

6 half drawer & 1 drawer unit




FI-A6A-30-7	30	762
FI-A6A-36-7	36	914
FI-A6A-42-7	42	1067
FI-A6A-48-7	48	1219

4 half drawer, 3 drawer & 1 door unit



FI-A7-15-8	15	381
FI-A7-18-8	18	457
FI-A7-21-8	21	533
FI-A7-24-8	24	610
FI-A7-30-8 *1	30	762
FI-A7-36-8 *1	36	914
FI-A7-42-8 *1*3	42	1067
FI-A7-48-8 *1*3	48	1219

8 half drawer unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FI-A8-30-5 *1	30	762
FI-A8-36-5 *1	36	914
FI-A8-42-5 *1*3	42	1067
FI-A8-48-5 *1*3	48	1219

5 drawer unit




Model no.	Width	
	inches	mm
FI-A10-30-10	30	762
FI-A10-36-10	36	914
FI-A10-42-10	42	1067
FI-A10-48-10	48	1219

4 half drawer & 6 drawer unit




FI-A8-30-7 *1	30	762
FI-A8-36-7 *1	36	914
FI-A8-42-7 *1*3	42	1067
FI-A8-48-7 *1*3	48	1219

4 half drawer & 3 drawer unit




FI-A10-30-12	30	762
FI-A10-36-12	36	914
FI-A10-42-12	42	1067
FI-A10-48-12	48	1219

8 half drawer & 4 drawer unit




FI-A9-30-6 *1	30	762
FI-A9-36-6 *1	36	914
FI-A9-42-6 *1*3	42	1067
FI-A9-48-6 *1*3	48	1219

6 drawer unit




FI-A10-30-14	30	762
FI-A10-36-14	36	914
FI-A10-42-14	42	1067
FI-A10-48-14	48	1219

12 half drawer & 2 drawer unit




FI-A9-30-10 *1	30	762
FI-A9-36-10 *1	36	914
FI-A9-42-10 *1*3	42	1067
FI-A9-48-10 *1*3	48	1219

8 half drawer & 2 drawer unit




FI-A10-30-16	30	762
FI-A10-36-16	36	914
FI-A10-42-16	42	1067
FI-A10-48-16	48	1219

16 half drawer unit




FI-A10-30-8	30	762
FI-A10-36-8	36	914
FI-A10-42-8	42	1067
FI-A10-48-8	48	1219

8 drawer unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FI-A11-15-2 *2	15	381
FI-A11-18-2 *2	18	457
FI-A11-21-2 *2	21	533
FI-A11-24-2 *2	24	610

2 file drawer unit




Model no.	Width	
	inches	mm
FI-A11-30-6 *2	30	762
FI-A11-36-6 *2	36	914
FI-A11-42-6 *2	42	1067
FI-A11-48-6 *2	48	1219

4 drawer & 2 file drawer unit




FI-A11-15-3 *2	15	381
FI-A11-18-3 *2	18	457
FI-A11-21-3 *2	21	533
FI-A11-24-3 *2	24	610

2 drawer & 1 file drawer unit




FI-A11-30-8 *2	30	762
FI-A11-36-8 *2	36	914
FI-A11-42-8 *2	42	1067
FI-A11-48-8 *2	48	1219

4 half drawer, 2 drawer & 2 file drawer unit




FI-A11-15-4 *2	15	381
FI-A11-18-4 *2	18	457
FI-A11-21-4 *2	21	533
FI-A11-24-4 *2	24	610

2 half drawer, 1 drawer & 1 file drawer unit




FI-A11-30-10 *2	30	762
FI-A11-36-10 *2	36	914
FI-A11-42-10 *2	42	1067
FI-A11-48-10 *2	48	1219

8 half drawer & 2 file drawer unit




FI-A11-15-5 *2	15	381
FI-A11-18-5 *2	18	457
FI-A11-21-5 *2	21	533
FI-A11-24-5 *2	24	610

4 half drawer & 1 file drawer unit




FI-A12-30-3 *2	30	762
FI-A12-36-3 *2	36	914
FI-A12-42-3 *2	42	1067
FI-A12-48-3 *2	48	1219

2 drawer, 1 file drawer & 1 door unit




FI-A11-30-4 *2	30	762
FI-A11-36-4 *2	36	914
FI-A11-42-4 *2	42	1067
FI-A11-48-4 *2	48	1219

4 file drawer unit



FI-A12-30-4 *2	30	762
FI-A12-36-4 *2	36	914
FI-A12-42-4 *2	42	1067
FI-A12-48-4 *2	48	1219

2 half drawer, 1 drawer, 1 file drawer & 1 door unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Standing Height Base Cabinets




Model no.	Width	
	inches	mm
FI-A12-30-5 *2	30	762
FI-A12-36-5 *2	36	914
FI-A12-42-5 *2	42	1067
FI-A12-48-5 *2	48	1219

4 half drawer, 1 file drawer & 1 door unit




FI-A13-30-4 *2	30	762
FI-A13-36-4 *2	36	914
FI-A13-42-4 *2	42	1067
FI-A13-48-4 *2	48	1219

3 drawer, 1 file drawer & 1 door unit




FI-A13-30-5 *2	30	762
FI-A13-36-5 *2	36	914
FI-A13-42-5 *2	42	1067
FI-A13-48-5 *2	48	1219

2 half drawer, 2 drawer, 1 file drawer & 1 door unit




FI-A13-30-6 *2	30	762
FI-A13-36-6 *2	36	914
FI-A13-42-6 *2	42	1067
FI-A13-48-6 *2	48	1219

4 half drawer, 1 drawer, 1 file drawer & 1 door unit



FI-A13-30-7 *2	30	762
FI-A13-36-7 *2	36	914
FI-A13-42-7 *2	42	1067
FI-A13-48-7 *2	48	1219

6 half drawer, 1 file drawer & 1 door unit




Model no.	Width	
	inches	mm
FI-A18-31	14,5	787

Description


- 2 pieces shelf.

45° corner unit




FI-A18S-31 *4	14,5	787
---------------	------	-----

45° corner sink unit




FI-A20-30	30	762
FI-A20-36	36	914
FI-A20-42	42	1067
FI-A20-48	48	1219

Corner unit



FI-A21-30	30	762
FI-A21-36	36	914
FI-A21-42	42	1067
FI-A21-48	48	1219


Corner unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.



Model no.	Width	
	inches	mm
FI-A15-15 *4	15	381
FI-A15-18 *4	18	457
FI-A15-21 *4	21	533
FI-A15-24 *4	24	610

1 door sink unit




FI-A15A-15 *4	15	381
FI-A15A-18 *4	18	457
FI-A15A-21 *4	21	533
FI-A15A-24 *4	24	610

Description


- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

1 door sink unit - false panel above



FI-A16-30 *4	30	762
FI-A16-36 *4	36	914
FI-A16-42 *4	42	1067
FI-A16-48 *4	48	1219

2 doors sink unit




FI-A16A-30 *4	30	762
FI-A16A-36 *4	36	914
FI-A16A-42 *4	42	1067
FI-A16A-48 *4	48	1219

Description


- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

2 door sink unit - false panel above



FI-A17-52 *4	52	1321
FI-A17-58 *4	58	1473

Island end sink unit




FI-A24-18	18	457
FI-A24-24	24	610

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

Remote control unit - 1 door




FI-A24-30	30	762
FI-A24-36	36	914
FI-A24-48	48	1219

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

Remote control unit - 2 doors




FI-A28-18	18	457
FI-A28-24	24	316

Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

1 door bin unit




FI-A29-36	36	914
FI-A29-48	48	1219

Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.


2 door bin unit



Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Standing Height Base Cabinets




Open unit

Model no.	Width	
	inches	mm
FI-A35-15	15	381
FI-A35-18	18	457
FI-A35-21	21	533
FI-A35-24	24	610
FI-A35-30	30	762
FI-A35-36	36	914
FI-A35-42	42	1067
FI-A35-48	48	1219

Description

- Two (2) adjustable shelves.




Pump unit, 1 door

FI-A41-15	15	381
FI-A41-18	18	457
FI-A41-21	21	533
FI-A41-24	24	610

Description


All pump cabinets are supplied with the following:

- Pull out tray located at the bottom of cabinet on 200 lbs capacity full extension tracks
- Prewired duplex 120V/20A plug inside cabinet with power switch located on front apron and wired to a junction box located at the back of the cabinet (access panel from inside the cabinet)
- Double wall perforated cabinet sides, back and top with sound insulation
- Unit can be supplied with "Silk screen" signage by adding "SC" after model number.




Pump unit, 2 doors

FI-A41-30	30	762
FI-A41-36	36	914
FI-A41-42	42	1067
FI-A41-48	48	1219




Pump Cabinet Signage

BLACK LETTERS




Door cabinet with stainless steel chute 8" x 6 1/2" high

FI-A42-15	15	381
FI-A42-18	18	457
FI-A42-21	21	533
FI-A42-24	24	610




Door cabinet with stainless steel chute 8" x 6 1/2" high

FI-A42-30	30	762
FI-A42-36	36	914
FI-A42-42	42	1067
FI-A42-48	48	1219




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Dimensions (inches)		
Model no.	width	depth
FI-A23-1821	18	21
FI-A23-2118	21	18
FI-A23-2121	21	21
FI-A23-2418	24	18
FI-A23-2421	24	21

Flammable unit - 1 door




Dimensions (inches)		
Model no.	width	depth
FI-A23-1821-CT	18	21
FI-A23-2118-CT	21	18
FI-A23-2121-CT	21	21
FI-A23-2418-CT	24	18
FI-A23-2421-CT	24	21

Mobile Flammable unit - 1 door



FI-A23-3018	30	18
FI-A23-3021	30	21
FI-A23-3618	36	18
FI-A23-3621	36	21
FI-A23-4218	42	18
FI-A23-4221	42	21
FI-A23-4818	48	18
FI-A23-4821	48	21

Flammable unit - 2 doors



FI-A23-3018-CT	30	18
FI-A23-3021-CT	30	21
FI-A23-3618-CT	36	18
FI-A23-3621-CT	36	21
FI-A23-4218-CT	42	18
FI-A23-4221-CT	42	21
FI-A23-4818-CT	48	18
FI-A23-4821-CT	48	21


Mobile Flammable unit - 2 doors

Description

For all flammable cabinets shown on this page:


- Solvent storage units meet O.S.H.A. standards and comply with N.F.P.A. code no 30 for flammable and combustible liquids
- An automatic door-closing mechanism activated in the event of fire is integrated and concealed within the cabinet structure. The doors may be held open by means of a chain equipped with a thermal fusible link.
- Two fire baffle vents are located in back of unit.
- Door is hinged with a continuous painted steel hinge
- Recessed base to match standard base cabinets.
- One galvanized steel shelf adjustable on 2" increment
- Painted steel liquid proof bottom pan 2" deep
- Three point locking handle
- Standard safety yellow exterior color
- Standard signage is applied with a sticker. Unit can be supplied with "Silk screen" by adding "SC" after model number
- Add "IP" at the end of the product code to order a removable liquid-tight painted steel tray to facilitate cleaning incase of spill.

Listed by Intertek to UL1275 and UL1275C




Flammable Cabinet Signage

RED LETTERS



Standing Height Corrosive Cabinets




Model no.	Width	
	inches	mm
FI-A25-18	18	457
FI-A25-24	24	610

Corrosive unit - 1 door

Description


For all corrosive cabinets shown on this page

- Unit is lined with chemical resistant molded polypropylene liner with an adjustable 2-positions HDPE shelf the full depth and width of the cabinet.
- Unit can be vented into fume hood by adding "V" after model number
- Unit can be supplied with "Silk screen" signage by adding "SC" after model number



FI-A25-18-18	18	457
FI-A25-24-18	24	610

Corrosive unit 18" deep - 1 door




FI-A25-30-SS	30	762
--------------	----	-----

Description

- Unit is supplied with four 4" high swivel casters with locks on the 2 front casters


*** Cabinets come with finished back and finished top**

Mobile Corrosive unit - 2 doors *




FI-A26-18	18	457
FI-A26-24	24	610

Corrosive unit - 1 door, false panel above




FI-A25-30	30	762
FI-A25-36	36	914
FI-A25-48	48	1219

Corrosive unit - 2 doors




FI-A26-18-18	18	457
FI-A26-24-18	24	610

Corrosive unit 18" deep - 1 door, false panel above



FI-A25-30-18	30	762
FI-A25-36-18	36	914
FI-A25-48-18	48	1219

Corrosive unit 18" deep - 2 doors




FI-A25-24-SS	24	610
--------------	----	-----

Description

- Unit is supplied with four 4" high swivel casters with locks on the 2 front casters


*** Cabinets come with finished back and finished top**

Mobile Corrosive unit - 1 door*




FI-A26-30	30	762
FI-A26-36	36	914
FI-A26-48	48	1219

Corrosive unit - 2 doors, blank panel above




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Corrosive Cabinet Signage

BLUE LETTERS




Model no.	Width	
	inches	mm
FI-B1-12-18	12	305
FI-B1-12-21	12	305
FI-B1-15	15	381
FI-B1-18	18	457
FI-B1-21	21	533
FI-B1-24	24	610

1 door unit




FI-B2-30	30	762
FI-B2-36	36	914
FI-B2-42	42	1067
FI-B2-48	48	1219

2 door unit



FI-B3-30	30	762
FI-B3-36	36	914
FI-B3-42	42	1067
FI-B3-48	48	1219

2 sliding door unit




FI-B4-15-1	15	381
FI-B4-18-1	18	457
FI-B4-21-1	21	533
FI-B4-24-1	24	610

1 drawer & 1 door unit




FI-B4-15-2	15	381
FI-B4-18-2	18	457
FI-B4-21-2	21	533
FI-B4-24-2	24	610

2 half drawer & 1 door unit




Model no.	Width	
	inches	mm
FI-B5-30-1 *1	30	762
FI-B5-36-1 *1	36	914
FI-B5-42-1 *1*3	42	1067
FI-B5-48-1 *1*3	48	1219

1 drawer & 2 door unit




FI-B5-30-2	30	762
FI-B5-36-2	36	914
FI-B5-42-2	42	1067
FI-B5-48-2	48	1219

2 drawer & 2 door unit




FI-B5-30-4	30	762
FI-B5-36-4	36	914
FI-B5-42-4	42	1067
FI-B5-48-4	48	1219

4 half drawer & 2 door unit




FI-B6-30-3	30	762
FI-B6-36-3	36	914
FI-B6-42-3	42	1067
FI-B6-48-3	48	1219

3 drawer & 1 door unit



FI-B6-30-4	30	762
FI-B6-36-4	36	914
FI-B6-42-4	42	1067
FI-B6-48-4	48	1219

2 half drawer, 2 drawer & 1 door unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

Sitting Height Base Cabinets




Model no.	Width	
	inches	mm
FI-B6-30-5	30	762
FI-B6-36-5	36	914
FI-B6-42-5	42	1067
FI-B6-48-5	48	1219

4 half drawer, 1 drawer & 1 door unit




FI-B6-30-6	30	762
FI-B6-36-6	36	914
FI-B6-42-6	42	1067
FI-B6-48-6	48	1219

6 half drawer unit




FI-B6A-30-4	30	762
FI-B6A-36-4	36	914
FI-B6A-42-4	42	1067
FI-B6A-48-4	48	1219

4 drawer & 1 door unit



FI-B7-15-3	15	381
FI-B7-18-3	18	457
FI-B7-21-3	21	533
FI-B7-24-3	24	610
FI-B7-30-3 *1	30	762
FI-B7-36-3 *1	36	914
FI-B7-42-3 *1*3	42	1067
FI-B7-48-3 *1*3	48	1219

3 drawer unit




FI-B7-15-4	15	381
FI-B7-18-4	18	457
FI-B7-21-4	21	533
FI-B7-24-4	24	610
FI-B7-30-4 *1	30	762
FI-B7-36-4 *1	36	914
FI-B7-42-4 *1*3	42	1067
FI-B7-48-4 *1*3	48	1219

2 half drawer & 2 drawer unit




Model no.	Width	
	inches	mm
FI-B7-15-6	15	381
FI-B7-18-6	18	457
FI-B7-21-6	21	533
FI-B7-24-6	24	610
FI-B7-30-6 *1	30	762
FI-B7-36-6 *1	36	914
FI-B7-42-6 *1*3	42	1067
FI-B7-48-6 *1*3	48	1219

6 half drawer unit




FI-B8-30-4 *1	30	762
FI-B8-36-4 *1	36	914
FI-B8-42-4 *1*3	42	1067
FI-B8-48-4 *1*3	48	1219

4 drawer unit




FI-B8-30-6 *1	30	762
FI-B8-36-6 *1	36	914
FI-B8-42-6 *1*3	42	1067
FI-B8-48-6 *1*3	48	1219

4 half drawer & 2 drawer unit




FI-B10-30-6	30	762
FI-B10-36-6	36	914
FI-B10-42-6	42	1067
FI-B10-48-6	48	1219

6 drawer unit



FI-B10-30-8	30	762
FI-B10-36-8	36	914
FI-B10-42-8	42	1067
FI-B10-48-8	48	1219

4 half drawer & 4 drawer unit



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Reverse models are also available by adding "RV" to model number.



Model no.	Width	
	inches	mm
FI-B10-30-12	30	762
FI-B10-36-12	36	914
FI-B10-42-12	42	1067
FI-B10-48-12	48	1219

12 half drawer unit




FI-B11-15-2 *2	15	381
FI-B11-18-2 *2	18	457
FI-B11-21-2 *2	21	533
FI-B11-24-2 *2	24	610

1 drawer & 1 file drawer unit




FI-B11-15-3 *2	15	381
FI-B11-18-3 *2	18	457
FI-B11-21-3 *2	21	533
FI-B11-24-3 *2	24	610

2 half drawer & 1 file drawer unit




FI-B11-30-4 *2	30	762
FI-B11-36-4 *2	36	914
FI-B11-42-4 *2	42	1067
FI-B11-48-4 *2	48	1219

2 drawer & 2 file drawer unit




FI-B11-30-6 *2	30	762
FI-B11-36-6 *2	36	914
FI-B11-42-6 *2	42	1067
FI-B11-48-6 *2	48	1219

4 half drawer & 2 file drawer unit




Model no.	Width	
	inches	mm
FI-B42-15	15	381
FI-B42-18	18	457
FI-B42-21	21	533
FI-B42-24	24	610

Door cabinet with stainless steel chute 8" x 6 1/2" w high



FI-B42-30	30	762
FI-B42-36	36	914
FI-B42-42	42	1067
FI-B42-48	48	1219


Door cabinet with stainless steel chute 8" x 6 1/2" high



Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

Sitting Height Base Cabinets




Model no.	Width	
	inches	mm
FI-B18-31	14,5	787

45° corner unit




FI-B18S-31 *4	14,5	787
---------------	------	-----

45° corner sink unit




FI-B20-30	30	762
FI-B20-36	36	914
FI-B20-42	42	1067
FI-B20-48	48	1219

Corner unit



FI-B21-30	30	762
FI-B21-36	36	914
FI-B21-42	42	1067
FI-B21-48	48	1219

Corner unit




FI-B28-18	18	457
FI-B28-24	24	316


Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

1 door bin unit




Open unit



Model no.	Width	
	inches	mm
FI-B35-15	15	381
FI-B35-18	18	457
FI-B35-21	21	533
FI-B35-24	24	610
FI-B35-30	30	762
FI-B35-36	36	914
FI-B35-42	42	1067
FI-B35-48	48	1219


FI-B15-15 *4	15	381
FI-B15-18 *4	18	457
FI-B15-21 *4	21	533
FI-B15-24 *4	24	610

1 door sink unit




FI-B16-30 *4	30	762
FI-B16-36 *4	36	914
FI-B16-42 *4	42	1067
FI-B16-48 *4	48	1219

2 door sink unit




FI-B17-52	52	1321
FI-B17-58	58	1473

Island end sink unit




Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

	Dimensions (inches)		
	Model no.	width	depth
	FI-B23-2418	24	18
	FI-B23-2421	24	21


Flammable unit - 1 door

	FI-B23-3018	30	18
	FI-B23-3021	30	21
	FI-B23-3618	36	18
	FI-B23-3621	36	21
	FI-B23-4218	42	18
	FI-B23-4221	42	21
	FI-B23-4818	48	18
	FI-B23-4821	48	21

Flammable unit - 2 doors

	Dimensions (inches)		
	Model no.	width	depth
	FI-B23-2418-CT	24	18
	FI-B23-2421-CT	24	21

Mobile Flammable unit - 1 door

	FI-B23-3018-CT	30	18
	FI-B23-3021-CT	30	21
	FI-B23-3618-CT	36	18
	FI-B23-3621-CT	36	21
	FI-B23-4218-CT	42	18
	FI-B23-4221-CT	42	21
	FI-B23-4818-CT	48	18
	FI-B23-4821-CT	48	21


Mobile Flammable unit - 2 doors

Description

For all flammable cabinets shown on this page:


- Solvent storage units meet O.S.H.A. standards and comply with N.F.P.A. code no 30 for flammable and combustible liquids
- An automatic door-closing mechanism activated in the event of fire is integrated and concealed within the cabinet structure. The doors may be held open by means of a chain equipped with a thermal fusible link.
- Two fire baffle vents are located in back of unit.
- Door is hinged with a continuous painted steel hinge
- Recessed base to match standard base cabinets.
- One galvanized steel shelf adjustable on 2" increment
- Painted steel liquid proof bottom pan 2" deep
- Three point locking handle
- Standard safety yellow exterior color
- Standard signage is applied with a sticker. Unit can be supplied with "Silk screen" by adding "SC" after model number
- Add "IP" at the end of the product code to order a removable liquid-tight painted steel tray to facilitate cleaning incase of spill.

Listed by Intertek to UL1275 and UL1275C




Flammable Cabinet Signage

RED LETTERS



Sitting Height Corrosive Cabinets




Model no.	Width	
	inches	mm
FI-B25-18	18	457
FI-B25-24	24	610

Description


For all corrosive cabinets shown on this page

- Unit is lined with a chemically resistant molded polypropylene liner with an adjustable 2-positions HDPE shelf the full depth and width of the cabinet.
- Unit can be vented into fume hood by adding "V" after model number
- Unit can be supplied with "Silk screen" signage by adding "SC" after model number

Corrosive unit - 1 door




FI-B25-18-18	18	457
FI-B25-24-18	24	610



FI-B25-30-18	30	762
FI-B25-36-18	36	914
FI-B25-48-18	48	1219

Corrosive unit 18" deep - 1 door



FI-B25-24-SS	24	533
--------------	----	-----


Description

- Unit is supplied with four swivel 4" high casters with locks on the 2 front casters

Cabinets come with finished back and finished top

Corrosive unit 18" deep - 2 doors

Mobile Corrosive unit - 1 door*




FI-B25-30-SS	30	610
--------------	----	-----

Description

- Unit is supplied with four swivel 4" high casters with locks on the 2 front casters


Cabinets come with finished back and finished top

Mobile Corrosive unit - 2 doors




FI-B25-30	30	762
FI-B25-36	36	914
FI-B25-48	48	1219

Corrosive unit - 2 doors




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Corrosive Cabinet Signage
BLUE LETTERS




Model no.	Width	
	inches	mm
FI-ADA1-12-18	12	305
FI-ADA1-12-21	12	305
FI-ADA1-15	15	381
FI-ADA1-18	18	457
FI-ADA1-21	21	533
FI-ADA1-24	24	610

1 door unit




FI-ADA2-30	30	762
FI-ADA2-36	36	914
FI-ADA2-42	42	1067
FI-ADA2-48	48	1219

2 door unit




FI-ADA3-30	30	762
FI-ADA3-36	36	914
FI-ADA3-42	42	1067
FI-ADA3-48	48	1219

2 sliding door unit




FI-ADA4-15-1	15	381
FI-ADA4-18-1	18	457
FI-ADA4-21-1	21	533
FI-ADA4-24-1	24	610

1 top drawer & 1 door unit




FI-ADA5-30-1	30	762
FI-ADA5-36-1	36	914
FI-ADA5-42-1 *3	42	1067
FI-ADA5-48-1 *3	48	1219

1 top drawer & 2 doors unit




Model no.	Width	
	inches	mm
FI-ADA5-30-2	30	762
FI-ADA5-36-2	36	914
FI-ADA5-42-2	42	1067
FI-ADA5-48-2	48	1219

2 top drawer & 2 door unit




FI-ADA6-30-4	30	762
FI-ADA6-36-4	36	914
FI-ADA6-42-4	42	1067
FI-ADA6-48-4	48	1219

3 drawer, 1 top drawer & 1 door unit




FI-ADA6A-30-5	30	762
FI-ADA6A-36-5	36	914
FI-ADA6A-42-5	42	1067
FI-ADA6A-48-5	48	1219

3 drawer, 2 top drawer & 1 door unit




FI-ADA6-30-7	30	762
FI-ADA6-36-7	36	914
FI-ADA6-42-7	42	1067
FI-ADA6-48-7	48	1219

6 drawer, 1 top drawer & 1 door unit



FI-ADA7-15-4	15	381
FI-ADA7-18-4	18	457
FI-ADA7-21-4	21	533
FI-ADA7-24-4	24	610
FI-ADA7-30-4 *1	30	762
FI-ADA7-36-4 *1	36	914
FI-ADA7-42-4 *1 *3	42	1067
FI-ADA7-48-4 *1 *3	48	1219

3 drawer & 1 top drawer unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high


Reverse models are also available by adding "RV" to model number.

ADA Height Base Cabinets




Model no.	Width	
	inches	mm
FI-ADA8-30-5	30	762
FI-ADA8-36-5	36	914
FI-ADA8-42-5 *3	42	1067
FI-ADA8-48-5 *3	48	1219

3 drawer & 2 top drawer unit




FI-ADA9-30-6	30	762
FI-ADA9-36-6	36	914
FI-ADA9-42-6 *3	42	1067
FI-ADA9-48-6 *3	48	1219

4 drawer & 2 top drawer unit




FI-ADA10-30-8	30	762
FI-ADA10-36-8	36	914
FI-ADA10-42-8	42	1067
FI-ADA10-48-8	48	1219

6 drawer & 2 top drawer unit




FI-ADA11-15-3 *2	15	381
FI-ADA11-18-3 *2	18	457
FI-ADA11-21-3 *2	21	533
FI-ADA11-24-3 *2	24	610

1 drawer, 1 top drawer & 1 file drawer unit




FI-ADA11-30-6 *2	30	762
FI-ADA11-36-6 *2	36	914
FI-ADA11-42-6 *2	42	1067
FI-ADA11-48-6 *2	48	1219

2 drawer, 2 top drawer & 2 file drawer unit



Model no.	Width	
	inches	mm
FI-ADA12-30-3 *2	30	762
FI-ADA12-36-3 *2	36	914
FI-ADA12-42-3 *2	42	1067
FI-ADA12-48-3 *2	48	1219


1 drawer, 1 top drawer, 1 file drawer & 1 door unit



Reverse models are also available by adding "RV" to model number.


Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Model no.	Width	
	inches	mm
FI-ADA18-31	14,5	787

45° corner unit




FI-ADA18S-31 *4	14,5	787
-----------------	------	-----

45° corner sink unit




FI-ADA20-30	30	762
FI-ADA20-36	36	914
FI-ADA20-42	42	1067
FI-ADA20-48	48	1219

Corner unit




FI-ADA21-30	30	762
FI-ADA21-36	36	914
FI-ADA21-42	42	1067
FI-ADA21-48	48	1219

Corner unit




Model no.	Width	
	inches	mm
FI-ADA15-15 *4	15	381
FI-ADA15-18 *4	18	457
FI-ADA15-21 *4	21	533
FI-ADA15-24 *4	24	610

1 door sink unit




FI-ADA15A-15 *4	15	381
FI-ADA15A-18 *4	18	457
FI-ADA15A-21 *4	21	533
FI-ADA15A-24 *4	24	610

1 door sink unit - false panel above




FI-ADA16-30 *4	30	762
FI-ADA16-36 *4	36	914
FI-ADA16-42 *4	42	1067
FI-ADA16-48 *4	48	1219

2 door sink unit




FI-ADA16A-30 *4	30	762
FI-ADA16A-36 *4	36	914
FI-ADA16A-42 *4	42	1067
FI-ADA16A-48 *4	48	1219

2 door sink unit - false panel above



FI-ADA17-52 *4	52	1321
FI-ADA17-58 *4	58	1473


Island end sink unit



Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

ADA Height Base Cabinets




Model no.	Width	
	inches	mm
FI-ADA24-18	18	457
FI-ADA24-24	24	610

Description

- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

Remote control unit - 1 door




FI-ADA24-30	30	762
FI-ADA24-36	36	914
FI-ADA24-48	48	1219

Description


- Fixed upper panel for mounting of services.
- This unit is mostly used under fume hoods where services are deck mounted.

Remote control unit - 2 doors




Open unit Description

- Two (2) adjustable shelves.



FI-ADA41-15	15	381
FI-ADA41-18	18	457
FI-ADA41-21	21	533
FI-ADA41-24	24	610

Pump unit, 1 door




FI-ADA41-30	30	762
FI-ADA41-36	36	914
FI-ADA41-42	42	1067
FI-ADA41-48	48	1219

Description

All pump cabinets are supplied with the following:

- Pull out tray located at the bottom of cabinet on 200 lbs capacity full extension tracks
- Prewired duplex 120V/20A plug inside cabinet with power switch located on front apron and wired to a junction box located at the rear of the cabinet (access panel from inside the cabinet)
- Double wall perforated interior cabinet sides, back and top with sound insulation
- Unit can be supplied with "Silk screen" signage by adding "SC" after model number

Pump unit, 2 doors




Model no.	Width	
	inches	mm
FI-ADA28-18	18	457
FI-ADA28-24	24	316

Description

- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

1 door bin unit




FI-ADA29-36	36	914
FI-ADA29-48	48	1219

Description


- Bins are not included.
- Bins can be supplied by adding the following letters after catalog number:
 - "GB" for galvanized steel;
 - "SB" for stainless steel T.304;
 - "LP" for leak proof.

2 doors bin unit




Pump Cabinet Signage

BLACK LETTERS




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Dimensions (inches)		
Model no.	width	depth
FI-ADA23-1821	18	21
FI-ADA23-2118	21	18
FI-ADA23-2121	21	21
FI-ADA23-2418	24	18
FI-ADA23-2421	24	21

Flammable unit - 1 door




Dimensions (inches)		
Model no.	width	depth
FI-ADA23-1821-CT	18	21
FI-ADA23-2118-CT	21	18
FI-ADA23-2121-CT	21	21
FI-ADA23-2418-CT	24	18
FI-ADA23-2421-CT	24	21

Mobile Flammable unit - 1 door



FI-ADA23-3018	30	18
FI-ADA23-3021	30	21
FI-ADA23-3618	36	18
FI-ADA23-3621	36	21
FI-ADA23-4218	42	18
FI-ADA23-4221	42	21
FI-ADA23-4818	48	18
FI-ADA23-4821	48	21

Flammable unit - 2 doors



FI-ADA23-3018-CT	30	18
FI-ADA23-3021-CT	30	21
FI-ADA23-3618-CT	36	18
FI-ADA23-3621-CT	36	21
FI-ADA23-4218-CT	42	18
FI-ADA23-4221-CT	42	21
FI-ADA23-4818-CT	48	18
FI-ADA23-4821-CT	48	21


Mobile Flammable unit - 2 doors

Description

For all flammable cabinets shown on this page:



- Solvent storage units meet O.S.H.A. standards and comply with N.F.P.A. code no 30 for flammable and combustible liquids
- An automatic door-closing mechanism activated in the event of fire is integrated and concealed within the cabinet structure. The doors may be held open by means of a chain equipped with a thermal fusible link.
- Two fire baffle vents are located in back of unit.
- Door is hinged with a continuous painted steel hinge
- Recessed base to match standard base cabinets.
- One galvanized steel shelf adjustable on 2" increment
- Painted steel liquid proof bottom pan 2" deep
- Three point locking handle
- Standard safety yellow exterior color
- Standard signage is applied with a sticker. Unit can be supplied with "Silk screen" by adding "SC" after model number
- Add "IP" at the end of the product code to order a removable liquid-tight painted steel tray to facilitate cleaning incase of spill.

Listed by Intertek to UL1275 and UL1275C




Flammable Cabinet Signage

RED LETTERS




ADA Height Corrosive Cabinets




Model no.	Width	
	inches	mm
FI-ADA25-18	18	457
FI-ADA25-24	24	610

Corrosive unit - 1 door




FI-ADA25-18-18	18	457
FI-ADA25-24-18	24	610

Corrosive unit 18" deep - 1 door




FI-ADA26-18	18	457
FI-ADA26-24	24	610

Corrosive unit - 1 door, false panel above



FI-ADA26-18-18	18	457
FI-ADA26-24-18	24	610

Corrosive unit 18" deep- 1 door, false panel above



FI-ADA25-24-SS	24	533
----------------	----	-----

Description

- Unit is supplied with four swivel 4" high casters with locks on the 2 front casters


Cabinets come with finished back and finished top

Mobile Corrosive unit - 1 door

Description

For all corrosive cabinets shown on this page

- Unit is lined with a chemically resistant molded polypropylene liner with an adjustable 2-positions HDPE shelf the full depth and width of the cabinet.
- Unit can be vented into fume hood by adding "V" after model number
- Unit can be supplied with "Silk screen" signage by adding "SC" after model number




FI-ADA25-30-SS	30	610
----------------	----	-----

Description

- Unit is supplied with four swivel 4" high casters with locks on the 2 front casters


Cabinets come with finished back and finished top

Mobile Corrosive unit - 2 doors




FI-ADA25-30	30	762
FI-ADA25-36	36	914
FI-ADA25-48	48	1219

Corrosive unit - 2 doors




FI-ADA25-30-18	30	762
FI-ADA25-36-18	36	914
FI-ADA25-48-18	48	1219

Corrosive unit 18" deep - 2 doors




FI-ADA26-30	30	762
FI-ADA26-36	36	914
FI-ADA26-48	48	1219

Corrosive unit - 2 doors, false panel above




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high




Corrosive Cabinet Signage

BLUE LETTERS




Model no.	Width	
	inches	mm
FI-OB1-15	15	381
FI-OB1-18	18	457
FI-OB1-21	21	533
FI-OB1-24	24	610

Suspended 1 Hinged Door unit




FI-OB2-30	30	762
FI-OB2-36	36	914
FI-OB2-42	42	1067
FI-OB2-48	48	1219

Suspended 2 Hinged Door unit




FI-OB4-15-1 *2	15	381
FI-OB4-18-1 *2	18	457
FI-OB4-21-1 *2	21	533
FI-OB4-24-1 *2	24	610

Suspended 1 Drawer & 1 Hinged Door unit



FI-OB5-30-2	30	762
FI-OB5-36-2	36	914
FI-OB5-42-2	42	1067
FI-OB5-48-2	48	1219

Suspended 2 Drawer & 2 Hinged Door unit




FI-OB7-15-3	15	381
FI-OB7-18-3	18	457
FI-OB7-21-3	21	533
FI-OB7-24-3	24	610
FI-OB7-30-3	30	762
FI-OB7-36-3	36	914
FI-OB7-42-3	42	1067
FI-OB7-48-3	48	1219

Suspended 3 Drawer unit




FI-OB7-15-4	15	381
FI-OB7-18-4	18	457
FI-OB7-21-4	21	533
FI-OB7-24-4	24	610
FI-OB7-30-4	30	762
FI-OB7-36-4	36	914
FI-OB7-42-4	42	1067
FI-OB7-48-4	48	1219

Suspended 2 Half Drawer & 2 Drawer unit




Model no.	Width	
	inches	mm
FI-OB10-30-6	30	762
FI-OB10-36-6	36	914
FI-OB10-42-6	42	1067
FI-OB10-48-6	48	1219

Suspended 6 Drawer unit




FI-OB10-30-8	30	762
FI-OB10-36-8	36	914
FI-OB10-42-8	42	1067
FI-OB10-48-8	48	1219

Suspended 4 Half Drawer & 4 Regular Drawer unit




FI-OB11-15-2 *2	15	381
FI-OB11-18-2 *2	18	457
FI-OB11-21-2 *2	21	533
FI-OB11-24-2 *2	24	610

Suspended 1 Drawer & 1 File Drawer unit




FI-OB11-15-3 *2	15	381
FI-OB11-18-3 *2	18	457
FI-OB11-21-3 *2	21	533
FI-OB11-24-3 *2	24	610

Suspended 2 Half Drawer & 1 File Drawer unit



FI-OB35-15	15	381
FI-OB35-18	18	457
FI-OB35-21	21	533
FI-OB35-24	24	610
FI-OB35-30	30	762
FI-OB35-36	36	914
FI-OB35-42	42	1067
FI-OB35-48	48	1219


Suspended Open unit




Legend:


- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high

19 3/8" Hight Suspended Cabinets


	Width		
	Model no.	inches	mm
	FI-OC1-15	15	381
	FI-OC1-18	18	457
	FI-OC1-21	21	533
	FI-OC1-24	24	610

	Width		
	Model no.	inches	mm
	FI-OC10-30-6	30	762
	FI-OC10-36-6	36	914
	FI-OC10-42-6	42	1067
	FI-OC10-48-6	48	1219


Suspended 1 Hinged Door unit

	FI-OC2-30	30	762
	FI-OC2-36	36	914
	FI-OC2-42	42	1067
	FI-OC2-48	48	1219


Suspended 2 Half Drawer & 4 Drawer unit

	FI-OC11-15-2 *2	15	381
	FI-OC11-18-2 *2	18	457
	FI-OC11-21-2 *2	21	533
	FI-OC11-24-2 *2	24	610


Suspended 2 Hinged Door unit

	FI-OC7-15-3	15	381
	FI-OC7-18-3	18	457
	FI-OC7-21-3	21	533
	FI-OC7-24-3	24	610


Suspended 1 Half Drawer + 1 File Drawer unit

	FI-OC11-30-4 *2	30	762
	FI-OC11-36-4 *2	36	914
	FI-OC11-42-4 *2	42	1067
	FI-OC11-48-4 *2	48	1219

Suspended 1 Half Drawer & 2 Drawer unit


	FI-OC7-15-4	15	381
	FI-OC7-18-4	18	457
	FI-OC7-21-4	21	533
	FI-OC7-24-4	24	610

Suspended 2 Half Drawer + 2 File Drawer unit

	FI-OC35-15	15	381
	FI-OC35-18	18	457
	FI-OC35-21	21	533
	FI-OC35-24	24	610
	FI-OC35-30	30	762
	FI-OC35-36	36	914
	FI-OC35-42	42	1067
	FI-OC35-48	48	1219


Suspended 3 Half Drawer & 1 Drawer unit

Suspended Open unit




Legend:

- *1 Double pulls on 30" to 48" cabinet sizes
- *2 File drawers supplied with hanger bars and full extension slide arms to fit legal or letter size files
- *3 Full extension drawer slide arms supplied as standard on these cabinets
- *4 Back panels in these cabinets are 9" high



Open shelf tall cabinet

Model no.	Width	
	inches	mm
FI-F10-84-15-14	15	381
FI-F10-84-18-14	18	457
FI-F10-84-21-14	21	533
FI-F10-84-24-14	24	610
FI-F10-84-30-14	30	762
FI-F10-84-36-14	36	914
FI-F10-84-42-14	42	1067
FI-F10-84-48-14	48	1219
FI-F10-84-15-18	15	381
FI-F10-84-18-18	18	457
FI-F10-84-21-18	21	533
FI-F10-84-24-18	24	610
FI-F10-84-30-18	30	762
FI-F10-84-36-18	36	914
FI-F10-84-42-18	42	1067
FI-F10-84-48-18	48	1219
FI-F10-84-15-22	15	381
FI-F10-84-18-22	18	457
FI-F10-84-21-22	21	533
FI-F10-84-24-22	24	610
FI-F10-84-30-22	30	762
FI-F10-84-36-22	36	914
FI-F10-84-42-22	42	1067
FI-F10-84-48-22	48	1219




Model no.	Width	
	inches	mm
FI-F20-84-30-14	30	762
FI-F20-84-36-14	36	914
FI-F20-84-42-14	42	1067
FI-F20-84-48-14	48	1219
FI-F20-84-30-18	30	762
FI-F20-84-36-18	36	914
FI-F20-84-42-18	42	1067
FI-F20-84-48-18	48	1219
FI-F20-84-30-22	30	762
FI-F20-84-36-22	36	914
FI-F20-84-42-22	42	1067
FI-F20-84-48-22	48	1219

2 hinged door tall cabinet




FI-F21-84-15-14	15	381
FI-F21-84-18-14	18	457
FI-F21-84-21-14	21	533
FI-F21-84-24-14	24	610
FI-F21-84-15-18	15	381
FI-F21-84-18-18	18	457
FI-F21-84-21-18	21	533
FI-F21-84-24-18	24	610
FI-F21-84-15-22	15	381
FI-F21-84-18-22	18	457
FI-F21-84-21-22	21	533
FI-F21-84-24-22	24	610

2 hinged door tall cabinet




FI-F20-84-15-14	15	381
FI-F20-84-18-14	18	457
FI-F20-84-21-14	21	533
FI-F20-84-24-14	24	610
FI-F20-84-15-18	15	381
FI-F20-84-18-18	18	457
FI-F20-84-21-18	21	533
FI-F20-84-24-18	24	610
FI-F20-84-15-22	15	381
FI-F20-84-18-22	18	457
FI-F20-84-21-22	21	533
FI-F20-84-24-22	24	610

1 hinged door tall cabinet




FI-F21-84-30-14	30	762
FI-F21-84-36-14	36	914
FI-F21-84-42-14	42	1067
FI-F21-84-48-14	48	1219
FI-F21-84-30-18	30	762
FI-F21-84-36-18	36	914
FI-F21-84-42-18	42	1067
FI-F21-84-48-18	48	1219
FI-F21-84-30-22	30	762
FI-F21-84-36-22	36	914
FI-F21-84-42-22	42	1067
FI-F21-84-48-22	48	1219

4 hinged door tall cabinet




Tall cabinet available depths of 14", 18" and 22"

84" High Tall Cabinets




Model no.	Width	
	inches	mm
FI-F30-84-15-14	15	381
FI-F30-84-18-14	18	457
FI-F30-84-21-14	21	533
FI-F30-84-24-14	24	610
FI-F30-84-15-18	15	381
FI-F30-84-18-18	18	457
FI-F30-84-21-18	21	533
FI-F30-84-24-18	24	610
FI-F30-84-15-22	15	381
FI-F30-84-18-22	18	457
FI-F30-84-21-22	21	533
FI-F30-84-24-22	24	610

1 glazed hinged door tall cabinet




Model no.	Width	
	inches	mm
FI-F31-84-15-14	15	381
FI-F31-84-18-14	18	457
FI-F31-84-21-14	21	533
FI-F31-84-24-14	24	610
FI-F31-84-15-18	15	381
FI-F31-84-18-18	18	457
FI-F31-84-21-18	21	533
FI-F31-84-24-18	24	610
FI-F31-84-15-22	15	381
FI-F31-84-18-22	18	457
FI-F31-84-21-22	21	533
FI-F31-84-24-22	24	610

2 glazed hinged door tall cabinet



FI-F30-84-30-14	30	762
FI-F30-84-36-14	36	914
FI-F30-84-42-14	42	1067
FI-F30-84-48-14	48	1219
FI-F30-84-30-18	30	762
FI-F30-84-36-18	36	914
FI-F30-84-42-18	42	1067
FI-F30-84-48-18	48	1219
FI-F30-84-30-22	30	762
FI-F30-84-36-22	36	914
FI-F30-84-42-22	42	1067
FI-F30-84-48-22	48	1219

2 glazed hinged door tall cabinet




FI-F31-84-30-14	30	762
FI-F31-84-36-14	36	914
FI-F31-84-42-14	42	1067
FI-F31-84-48-14	48	1219
FI-F31-84-30-18	30	762
FI-F31-84-36-18	36	914
FI-F31-84-42-18	42	1067
FI-F31-84-48-18	48	1219
FI-F31-84-30-22	30	762
FI-F31-84-36-22	36	914
FI-F31-84-42-22	42	1067
FI-F31-84-48-22	48	1219


4 glazed hinged door tall cabinet

Description

- Framed glass doors are supplied using 3 mm float glass




Tall cabinet: available depths of 14", 18" and 22"




Model no.	Width	
	inches	mm
FI-F40-84-30-14	30	762
FI-F40-84-36-14	36	914
FI-F40-84-42-14	42	1067
FI-F40-84-48-14	48	1219
FI-F40-84-30-18	30	762
FI-F40-84-36-18	36	914
FI-F40-84-42-18	42	1067
FI-F40-84-48-18	48	1219
FI-F40-84-30-22	30	762
FI-F40-84-36-22	36	914
FI-F40-84-42-22	42	1067
FI-F40-84-48-22	48	1219

2 sliding door tall cabinet




Model no.	Width	
	inches	mm
FI-F90-84-30-14	30	762
FI-F90-84-36-14	36	914
FI-F90-84-42-14	42	1067
FI-F90-84-48-14	48	1219
FI-F90-84-30-18	30	762
FI-F90-84-36-18	36	914
FI-F90-84-42-18	42	1067
FI-F90-84-48-18	48	1219
FI-F90-84-30-22	30	762
FI-F90-84-36-22	36	914
FI-F90-84-42-22	42	1067
FI-F90-84-48-22	48	1219

4 glass sliding door tall cabinet



FI-F41-84-30-14	30	762
FI-F41-84-36-14	36	914
FI-F41-84-42-14	42	1067
FI-F41-84-48-14	48	1219
FI-F41-84-30-18	30	762
FI-F41-84-36-18	36	914
FI-F41-84-42-18	42	1067
FI-F41-84-48-18	48	1219
FI-F41-84-30-22	30	762
FI-F41-84-36-22	36	914
FI-F41-84-42-22	42	1067
FI-F41-84-48-22	48	1219

4 sliding door tall cabinet




FI-F50-84-30-14	30	762
FI-F50-84-36-14	36	914
FI-F50-84-42-14	42	1067
FI-F50-84-48-14	48	1219
FI-F50-84-30-18	30	762
FI-F50-84-36-18	36	914
FI-F50-84-42-18	42	1067
FI-F50-84-48-18	48	1219
FI-F50-84-30-22	30	762
FI-F50-84-36-22	36	914
FI-F50-84-42-22	42	1067
FI-F50-84-48-22	48	1219

2 glazed sliding door tall cabinet


Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass



FI-F51-84-30-14	30	762
FI-F51-84-36-14	36	914
FI-F51-84-42-14	42	1067
FI-F51-84-48-14	48	1219
FI-F51-84-30-18	30	762
FI-F51-84-36-18	36	914
FI-F51-84-42-18	42	1067
FI-F51-84-48-18	48	1219
FI-F51-84-30-22	30	762
FI-F51-84-36-22	36	914
FI-F51-84-42-22	42	1067
FI-F51-84-48-22	48	1219

4 glazed sliding door tall cabinet




Tall cabinet available depths of 14", 18" and 22"

24" High Wall Cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FE10-24-15-14	15	381
FE10-24-18-14	18	457
FE10-24-21-14	21	533
FE10-24-24-14	24	610
FE10-24-30-14	30	762
FE10-24-36-14	36	914
FE10-24-42-14	42	1067
FE10-24-48-14	48	1219
FE10-24-15-18	15	381
FE10-24-18-18	18	457
FE10-24-21-18	21	533
FE10-24-24-18	24	610
FE10-24-30-18	30	762
FE10-24-36-18	36	914
FE10-24-42-18	42	1067
FE10-24-48-18	48	1219




1 hinged door unit

FE20-24-15-14	15	381
FE20-24-18-14	18	457
FE20-24-21-14	21	533
FE20-24-24-14	24	610
FE20-24-15-18	15	381
FE20-24-18-18	18	457
FE20-24-21-18	21	533
FE20-24-24-18	24	610




1 glazed hinged door unit

FE30-24-15-14	15	381
FE30-24-18-14	18	457
FE30-24-21-14	21	533
FE30-24-24-14	24	610
FE30-24-15-18	15	381
FE30-24-18-18	18	457
FE30-24-21-18	21	533
FE30-24-24-18	24	610




2 hinged door unit

FE20-24-30-14	30	762
FE20-24-36-14	36	914
FE20-24-42-14	42	1067
FE20-24-48-14	48	1219
FE20-24-30-18	30	762
FE20-24-36-18	36	914
FE20-24-42-18	42	1067
FE20-24-48-18	48	1219




2 glazed hinged door unit

Model no.	Width	
	inches	mm
FE30-24-30-14	30	762
FE30-24-36-14	36	914
FE30-24-42-14	42	1067
FE30-24-48-14	48	1219
FE30-24-30-18	30	762
FE30-24-36-18	36	914
FE30-24-42-18	42	1067
FE30-24-48-18	48	1219




2 glass sliding door unit

FE90-24-30-14	30	762
FE90-24-36-14	36	914
FE90-24-42-14	42	1067
FE90-24-48-14	48	1219
FE90-24-30-18	30	762
FE90-24-36-18	36	914
FE90-24-42-18	42	1067
FE90-24-48-18	48	1219



2 solid sliding door unit

FE40-24-30-14	30	762
FE40-24-36-14	36	914
FE40-24-42-14	42	1067
FE40-24-48-14	48	1219
FE40-24-30-18	30	762
FE40-24-36-18	36	914
FE40-24-42-18	42	1067
FE40-24-48-18	48	1219




2 glazed sliding door unit


FE50-24-30-14	30	762
FE50-24-36-14	36	914
FE50-24-42-14	42	1067
FE50-24-48-14	48	1219
FE50-24-30-18	30	762
FE50-24-36-18	36	914
FE50-24-42-18	42	1067
FE50-24-48-18	48	1219

Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




14" depth units




18" depth units

Reverse models are also available by adding "RV" to model number.




Model no.	Width	
	inches	mm
FI-E22-24-30-14	30	762
FI-E22-24-36-14	36	914
FI-E22-24-42-14	42	1067
FI-E22-24-48-14	48	1219
FI-E22-24-30-18	30	762
FI-E22-24-36-18	36	914
FI-E22-24-42-18	42	1067
FI-E22-24-48-18	48	1219

Corner unit with door




14" depth units 18" depth units

Reverse models are also available by adding "RV" to model number.




Model no.	Width	
	inches	mm
FI-E32-24-30-14	30	762
FI-E32-24-36-14	36	914
FI-E32-24-42-14	42	1067
FI-E32-24-48-14	48	1219
FI-E32-24-30-18	30	762
FI-E32-24-36-18	36	914
FI-E32-24-42-18	42	1067
FI-E32-24-48-18	48	1219

Corner unit with glazed door




14" depth units 18" depth units

Reverse models are also available by adding "RV" to model number.




FI-E23-24-14-14	14	365
FI-E23-24-18-18	18	457

45° corner unit with door



FI-E33-24-14-14	14	365
FI-E33-24-18-18	18	457

45° corner unit with glazed door




14" depth units 18" depth units

Description


- Framed glass doors are supplied using 3 mm float glass

30" High Wall Cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FE10-30-15-14	15	381
FE10-30-18-14	18	457
FE10-30-21-14	21	533
FE10-30-24-14	24	610
FE10-30-30-14	30	762
FE10-30-36-14	36	914
FE10-30-42-14	42	1067
FE10-30-48-14	48	1219
FE10-30-15-18	15	381
FE10-30-18-18	18	457
FE10-30-21-18	21	533
FE10-30-24-18	24	610
FE10-30-30-18	30	762
FE10-30-36-18	36	914
FE10-30-42-18	42	1067
FE10-30-48-18	48	1219




2 glazed hinged door unit

Model no.	Width	
	inches	mm
FE30-30-30-14	30	762
FE30-30-36-14	36	914
FE30-30-42-14	42	1067
FE30-30-48-14	48	1219
FE30-30-30-18	30	762
FE30-30-36-18	36	914
FE30-30-42-18	42	1067
FE30-30-48-18	48	1219




2 glass sliding door unit

FE90-30-30-14	30	762
FE90-30-36-14	36	914
FE90-30-42-14	42	1067
FE90-30-48-14	48	1219
FE90-30-30-18	30	762
FE90-30-36-18	36	914
FE90-30-42-18	42	1067
FE90-30-48-18	48	1219




1 hinged door unit

FE20-30-15-14	15	381
FE20-30-18-14	18	457
FE20-30-21-14	21	533
FE20-30-24-14	24	610
FE20-30-15-18	15	381
FE20-30-18-18	18	457
FE20-30-21-18	21	533
FE20-30-24-18	24	610




2 solid sliding door unit

FE40-30-30-14	30	762
FE40-30-36-14	36	914
FE40-30-42-14	42	1067
FE40-30-48-14	48	1219
FE40-30-30-18	30	762
FE40-30-36-18	36	914
FE40-30-42-18	42	1067
FE40-30-48-18	48	1219




1 glazed hinged door unit

FE30-30-15-14	15	381
FE30-30-18-14	18	457
FE30-30-21-14	21	533
FE30-30-24-14	24	610
FE30-30-15-18	15	381
FE30-30-18-18	18	457
FE30-30-21-18	21	533
FE30-30-24-18	24	610



2 glazed sliding door unit

FE50-30-30-14	30	762
FE50-30-36-14	36	914
FE50-30-42-14	42	1067
FE50-30-48-14	48	1219
FE50-30-30-18	30	762
FE50-30-36-18	36	914
FE50-30-42-18	42	1067
FE50-30-48-18	48	1219




2 hinged door unit


FE20-30-30-14	30	762
FE20-30-36-14	36	914
FE20-30-42-14	42	1067
FE20-30-48-14	48	1219
FE20-30-30-18	30	762
FE20-30-36-18	36	914
FE20-30-42-18	42	1067
FE20-30-48-18	48	1219

Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




14" depth units




18" depth units

Reverse models are also available by adding "RV" to model number.




Corner unit with door

Model no.	Width	
	inches	mm
FE22-30-30-14	30	762
FE22-30-36-14	36	914
FE22-30-42-14	42	1067
FE22-30-48-14	48	1219
FE22-30-30-18	30	762
FE22-30-36-18	36	914
FE22-30-42-18	42	1067
FE22-30-48-18	48	1219




Reverse models are also available by adding "RV" to model number.




Corner unit with glazed door

Model no.	Width	
	inches	mm
FE32-30-30-14	30	762
FE32-30-36-14	36	914
FE32-30-42-14	42	1067
FE32-30-48-14	48	1219
FE32-30-30-18	30	762
FE32-30-36-18	36	914
FE32-30-42-18	42	1067
FE32-30-48-18	48	1219




Reverse models are also available by adding "RV" to model number.




45° corner unit with door

FE23-30-14-14	14	365
FE23-30-18-18	18	457



45° corner unit with glazed door


FE33-30-14-14	14	365
FE33-30-18-18	18	457



Description


- Framed glass doors are supplied using 3 mm float glass

36" High Wall Cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FE10-36-15-14	15	381
FE10-36-18-14	18	457
FE10-36-21-14	21	533
FE10-36-24-14	24	610
FE10-36-30-14	30	762
FE10-36-36-14	36	914
FE10-36-42-14	42	1067
FE10-36-48-14	48	1219
FE10-36-15-18	15	381
FE10-36-18-18	18	457
FE10-36-21-18	21	533
FE10-36-24-18	24	610
FE10-36-30-18	30	762
FE10-36-36-18	36	914
FE10-36-42-18	42	1067
FE10-36-48-18	48	1219




2 glazed hinged door unit

Model no.	Width	
	inches	mm
FE30-36-30-14	30	762
FE30-36-36-14	36	914
FE30-36-42-14	42	1067
FE30-36-48-14	48	1219
FE30-36-30-18	30	762
FE30-36-36-18	36	914
FE30-36-42-18	42	1067
FE30-36-48-18	48	1219




2 glass sliding door unit

FE90-36-30-14	30	762
FE90-36-36-14	36	914
FE90-36-42-14	42	1067
FE90-36-48-14	48	1219
FE90-36-30-18	30	762
FE90-36-36-18	36	914
FE90-36-42-18	42	1067
FE90-36-48-18	48	1219




1 hinged door unit

FE20-36-15-14	15	381
FE20-36-18-14	18	457
FE20-36-21-14	21	533
FE20-36-24-14	24	610
FE20-36-15-18	15	381
FE20-36-18-18	18	457
FE20-36-21-18	21	533
FE20-36-24-18	24	610




2 sliding door unit

FE40-36-30-14	30	762
FE40-36-36-14	36	914
FE40-36-42-14	42	1067
FE40-36-48-14	48	1219
FE40-36-30-18	30	762
FE40-36-36-18	36	914
FE40-36-42-18	42	1067
FE40-36-48-18	48	1219




1 glazed hinged door unit

FE30-36-15-14	15	381
FE30-36-18-14	18	457
FE30-36-21-14	21	533
FE30-36-24-14	24	610
FE30-36-15-18	15	381
FE30-36-18-18	18	457
FE30-36-21-18	21	533
FE30-36-24-18	24	610



2 glazed sliding door unit

FE50-36-30-14	30	762
FE50-36-36-14	36	914
FE50-36-42-14	42	1067
FE50-36-48-14	48	1219
FE50-36-30-18	30	762
FE50-36-36-18	36	914
FE50-36-42-18	42	1067
FE50-36-48-18	48	1219




2 solid hinged door unit


FE20-36-30-14	30	762
FE20-36-36-14	36	914
FE20-36-42-14	42	1067
FE20-36-48-14	48	1219
FE20-36-30-18	30	762
FE20-36-36-18	36	914
FE20-36-42-18	42	1067
FE20-36-48-18	48	1219

Description

- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass




14" depth units




18" depth units

Reverse models are also available by adding "RV" to model number.




Corner unit with door




Reverse models are also available by adding "RV" to model number.

Model no.	Width	
	inches	mm
FE22-36-30-14	30	762
FE22-36-36-14	36	914
FE22-36-42-14	42	1067
FE22-36-48-14	48	1219
FE22-36-30-18	30	762
FE22-36-36-18	36	914
FE22-36-42-18	42	1067
FE22-36-48-18	48	1219




Corner unit with glazed door




Reverse models are also available by adding "RV" to model number.

Model no.	Width	
	inches	mm
FE32-36-30-14	30	762
FE32-36-36-14	36	914
FE32-36-42-14	42	1067
FE32-36-48-14	48	1219
FE32-36-30-18	30	762
FE32-36-36-18	36	914
FE32-36-42-18	42	1067
FE32-36-48-18	48	1219




45° corner unit with door

FE23-36-14-14	14	365
FE23-36-18-18	18	457




45° corner unit with glazed door

FE33-36-14-14	14	365
FE33-36-18-18	18	457




- Description**
- Framed glass doors are supplied using 3 mm float glass

48" High Wall Cabinets




Open shelf unit

Model no.	Width	
	inches	mm
FE10-48-15-14	15	381
FE10-48-18-14	18	457
FE10-48-21-14	21	533
FE10-48-24-14	24	610
FE10-48-30-14	30	762
FE10-48-36-14	36	914
FE10-48-42-14	42	1067
FE10-48-48-14	48	1219
FE10-48-15-18	15	381
FE10-48-18-18	18	457
FE10-48-21-18	21	533
FE10-48-24-18	24	610
FE10-48-30-18	30	762
FE10-48-36-18	36	914
FE10-48-42-18	42	1067
FE10-48-48-18	48	1219




2 glazed hinged door unit

Model no.	Width	
	inches	mm
FE30-48-30-14	30	762
FE30-48-36-14	36	914
FE30-48-42-14	42	1067
FE30-48-48-14	48	1219
FE30-48-30-18	30	762
FE30-48-36-18	36	914
FE30-48-42-18	42	1067
FE30-48-48-18	48	1219




2 glass sliding door unit

FE90-48-30-14	30	762
FE90-48-36-14	36	914
FE90-48-42-14	42	1067
FE90-48-48-14	48	1219
FE90-48-30-18	30	762
FE90-48-36-18	36	914
FE90-48-42-18	42	1067
FE90-48-48-18	48	1219




1 glazed hinged door unit

FE20-48-15-14	15	381
FE20-48-18-14	18	457
FE20-48-21-14	21	533
FE20-48-24-14	24	610
FE20-48-15-18	15	381
FE20-48-18-18	18	457
FE20-48-21-18	21	533
FE20-48-24-18	24	610




2 sliding door unit

FE40-48-30-14	30	762
FE40-48-36-14	36	914
FE40-48-42-14	42	1067
FE40-48-48-14	48	1219
FE40-48-30-18	30	762
FE40-48-36-18	36	914
FE40-48-42-18	42	1067
FE40-48-48-18	48	1219




1 hinged door unit

FE30-48-15-14	15	381
FE30-48-18-14	18	457
FE30-48-21-14	21	533
FE30-48-24-14	24	610
FE30-48-15-18	15	381
FE30-48-18-18	18	457
FE30-48-21-18	21	533
FE30-48-24-18	24	610



2 glazed sliding door unit

FE50-48-30-14	30	762
FE50-48-36-14	36	914
FE50-48-42-14	42	1067
FE50-48-48-14	48	1219
FE50-48-30-18	30	762
FE50-48-36-18	36	914
FE50-48-42-18	42	1067
FE50-48-48-18	48	1219




2 hinged door unit


FE20-48-30-14	30	762
FE20-48-36-14	36	914
FE20-48-42-14	42	1067
FE20-48-48-14	48	1219
FE20-48-30-18	30	762
FE20-48-36-18	36	914
FE20-48-42-18	42	1067
FE20-48-48-18	48	1219

Description


- Unframed glass sliding doors are supplied using 5 mm float glass with finger grooves
- Framed glass doors are supplied using 3 mm float glass



14" depth units




18" depth units




Corner unit with door

Model no.	Width	
	inches	mm
FlE22-48-30-14	30	762
FlE22-48-36-14	36	914
FlE22-48-42-14	42	1067
FlE22-48-48-14	48	1219
FlE22-48-30-18	30	762
FlE22-48-36-18	36	914
FlE22-48-42-18	42	1067
FlE22-48-48-18	48	1219




14" depth units 18" depth units

Reverse models are also available by adding "RV" to model number.




Corner unit with glazed door

Model no.	Width	
	inches	mm
FlE32-48-30-14	30	762
FlE32-48-36-14	36	914
FlE32-48-42-14	42	1067
FlE32-48-48-14	48	1219
FlE32-48-30-18	30	762
FlE32-48-36-18	36	914
FlE32-48-42-18	42	1067
FlE32-48-48-18	48	1219




14" depth units

Reverse models are also available by adding "RV" to model number.




45° corner unit with door

FlE23-48-14-14	14	365
FlE23-48-18-18	18	457



45° corner unit with glazed door

FlE33-48-14-14	14	365
FlE33-48-18-18	18	457




14" depth units

18" depth units

Description

- Framed glass doors are supplied using 3 mm float glass

1 CABINET PULLS**Standard — Brushed nickel wire pull 4" c/c**WP ☐ Stainless steel wire pull 4" c/cPP ☐ Black recessed 4" x 2" pullAP ☐ Aluminum recessed 4" x 2" pull**2 DRAWER SLIDE ARMS****Standard****75% opening self closing slide arms on all regular drawers – SEFA 8M – 150 lbs****100% full extension opening on file drawers – SEFA 8M – 200 lbs****100% full extension opening on all drawers from 36" to 48" wide**TZ ☐ 100% full extension opening on all project drawers – SEFA 8M – 100 lbsTY ☐ 100% full extension opening with self closing option on all project drawers – SEFA 8M – 100 lbsTX ☐ 100% full extension opening on all project drawers – SEFA 8M – 200 lbs**3 CABINET HINGES****Standard — 5-knuckles stainless steel vibra finished**BH ☐ 5-knuckles black finishSH ☐ 5-knuckles stainless steel**4 LOCKS****Standard — None**KA ☐ 8-disc chrome plated with removable core on all doors and all drawersKT ☐ 8-disc chrome plated with removable core on all drawersKP ☐ 8-disc chrome plated with removable core on all doorsK1 ☐ 8-disc chrome plated with removable core on top drawers onlyKF ☐ 8-disc chrome plated with removable core on file drawers onlyKC ☐ Chrome plated padlock bracket on all doors and all drawersXI ☐ Interlocking system that permits only one drawer opening at the time (with full extension arms only) – standard on mobile cabinetsTP ☐ Theft security panels between drawersL1 ☐ All locks different with a project master keyL2 ☐ All locks the same per cabinet with a project master keyL3 ☐ All locks the same per room with a room master keyL4 ☐ All locks the same**5 DOOR CATCHES****Standard — Spring loaded roller door catch****6 CARD HOLDERS****Standard — None**CH ☐ Aluminum 4" x 2" card holder on all drawersCD ☐ Aluminum 4" x 2" card holder on all doors

7 CABINET SHELF

Standard — One (1) full depth shelf, adjustable on ½" increment

- T1 ☐ Additional shelf per cabinet
- P1 ☐ Pull out shelf for cabinets 15", 18", 21" or 24" only
- PR ☐ Pull out shelf on right cabinet side
- PL ☐ Pull out shelf on left cabinet side
- S2 ☐ Split shelves (2 pieces), adjustable on ½" increment
- FS ☐ Cabinet full extension shelf supplied with full extension slide arms – 200 lbs

8 FINISHED PANELS

Standard — Removable access panel on cupboard cabinets only (open back on drawer cab)

- BP ☐ Removable access panel on drawer units
- FO ☐ Painted steel flat finished back welded to cabinet
- FT ☐ Painted steel flat finished top and back panels welded to cabinet

9 CABINET BASE OPTIONS

Standard — 3" x 4" high integral front toe kick with four 3/8"x 2 ½" leveler glides

- PN ☐ Cabinet levelers supplied with grey rubber cap
- SS ☐ Four swivel casters 4" height, rubber, with brakes on two front casters – 275 lbs
- WB ☐ Cabinet supplied with no toe kick (4" shorter)
- RT ☐ Cabinet supplied with a right toe kick 3" x 4" high – not available on tall cabinets
- LT ☐ Cabinet supplied with a left toe kick 3" x 4" high – not available on tall cabinets
- VT ☐ Cabinet supplied with a venting grill on the front toe kick

10 VARIOUS OTHER CABINET OPTIONS

- RV ☐ Reverse door/drawer arrangement
- SK ☐ Sink support consisting of 2 horizontal struts with 4 rod hangers
- LG ☐ 6mm Laminated glass in lieu of standard 3mm float glass
- TG ☐ 6mm Tempered glass in lieu of standard 3mm float glass
- SB ☐ Stainless steel bin insert into bin cabinets
- GB ☐ Galvanized bin insert into bin cabinets
- LP ☐ Leak proof bin insert into bin cabinets
- WC ☐ Louvers incorporated into cabinet door(s)
- SS4 ☐ Same cabinet but supplied in stainless steel T.304
- SS6 ☐ Same cabinet but supplied in stainless steel T.316



BEDCOLAB LTD

**2305, Francis-Hugues Avenue
Laval, Quebec, Canada H7S 1N5**

**Phone: 514.384.2820 / 1.800.461.6414
Fax: 514.384.9210**

**E-mail: information@bedcolab.com
Web Site: www.bedcolab.com**



Authorized dealer: