

(bedcolab)



MEZZO

Multi-purpose table

Built to perform. Anywhere.

More than just a work surface, this multi-purpose table is built to support ideas in motion. Whether it's shaping young minds in a K-12 lab or advancing innovation in a research facility, its sturdy and adaptable design provides the strength and flexibility every workspace deserves.



DATA SHEET
MEZZO



Dimensions (in)
Includes 1" countertop

Width: 30, 36, 42, 48, 54, 60, 66, 72
Depth: 24, 30, 36
Fixed Height: 36
Adjustable Height: 30 to 37 (1" increments)

Construction

All-welded structure
Legs, 2"x2" square tubing
Side and rear aprons, 4" high
Suspension rails, 1¼" high
Rear and side cross-member, 1"x2"
Levelers
Painted steel, 20 SEFA compliant colors

Options

Casters (four casters, front legs with brakes)
Stainless Steel 304

Countertop


Epoxy
Phenolic
Laminate

Certification

SEFA 8M (600 lbs max load capacity)

Related products

Forte Inset, Forte Flush Overlay and Harmony suspended cabinets. Aria Bench, both have square tubing legs



PRODUCT CONFIGURATIONS



ME-FX

Fixed height table with levelers



ME-AJ

Adjustable height table with levelers



ME-FX with CT option

Fixed height table on casters




ME-AJ with CT option

Adjustable height table on casters

PRODUCT CODE STRUCTURE

This page was created to explain how product codes are built and what is the meaning of every element in them. Using a standard product code in your project documents allows us to process your order quickly and accurately — so you get exactly what you expect, every time.

The complete list of product codes can be found in our [online library](#) under the “Product code lists” section.



Leg type

BB FX: Fixed Height
AJ: Adjustable Height

Dimensions (inches)

CC Width
DD Depth

Options

EE CT: Casters
FFF SS4: Stainless Steel 304

**Find a Bedcolab
dealer near you**

