

(bedcolab)

Laboratory Casework Systems and Fume Hoods



MEZZO

Multi-purpose table

Built to perform. Anywhere.

More than just a work surface, these multi-purpose tables are built to support ideas in motion. Whether it's shaping young minds in a K-12 lab or advancing innovation in a research facility, their sturdy and adaptable design provide the strength and flexibility every workspace deserves.



Emory University HSRB II
Executed in collaboration with H2I Group

DATA SHEET

MEZZO



Dimensions (in)
Includes 1" countertop

Width: 30, 36, 42, 48, 54, 60, 66, 72
Depth: 24, 30, 36
Fixed Height: 36
Adjustable Height: 30 to 37 (1" increments)

Construction

All-welded structure
Legs, 2"x2" square tubing
Side and rear aprons, 4" high
Cabinet suspension rails, 1¼" high
Rear and side cross-member, 1"x2"
Levelers
Painted steel, 20 SEFA compliant colors

Options

Casters (four casters, front legs with brakes)
Stainless Steel 304

Countertop
Provided by Dealer

Epoxy
Phenolic
Laminate


Certification

SEFA 8M (600 lbs max load capacity)

Related products

Forte Inset, Forte Flush Overlay and Harmony Flush Overlay suspended cabinets
• OC Series compatible with all Mezzo tables
• OB Series compatible with 30" and 36" deep Mezzo tables

Aria Bench and Mezzo Table have the same leg style



CONFIGURATIONS

MEZZO



ME-FX
Fixed height table with levelers



ME-AJ
Adjustable height table with levelers



ME-FX with CT option
Fixed height table on casters




ME-AJ with CT option
Adjustable height table on casters

PRODUCT CODE STRUCTURE

MEZZO

This page was created to explain how product codes are built and what is the meaning of every element in them. Using a standard product code in your project documents allows us to process your order quickly and accurately — so you get exactly what you expect, every time.

The complete list of product codes can be found in our [online library](#) under the “Product code lists” section.



| | |
|-------------------------------|--|
| Product line | AA ME: Mezzo |
| Leg type | BB FX: Fixed Height AJ: Adjustable Height |
| Dimensions <i>(inches)</i> | CC Width DD Depth |
| Options | EE CT: Casters FFF SS4: Stainless Steel 304 |

**Find a Bedcolab
dealer near you**

